

untitled protocol

Sunamita Lima

Abstract

For isolated diaphragm skeletal muscle preparation, each hemi-diaphragm segment was transferred to an organ bath, and the dissected tissue was tied up to a holder fixed at the bottom of the glass organ bath, while the upper extremity of the tissue was fixed to a force transducer. The contraction of the diaphragm muscle was induced by transmural electrical stimuli conducted through the platinum electrodes.

Citation: Sunamita Lima untitled protocol. protocols.io

dx.doi.org/10.17504/protocols.io.p7vdrn6

Published: 17 May 2018

Materials

- Sodium phosphate dibasic [7558-79-4](#) by [Sigma Aldrich](#)
- ✓ NaHCO₃ by Contributed by users
- ✓ KCl by Contributed by users
- ✓ CaCl₂ by Contributed by users
- NaCl [53014](#) by [Sigma Aldrich](#)
- ✓ MgCl₆.H₂O by Contributed by users
- ✓ d-glucose by Contributed by users

Protocol

1

Step 1.

The animal was euthanized by cervical dislocation

1

Step 2.

the diaphragm muscle was removed and transferred to a petri dish containing Tyrode's solution

1

Step 3.

A 1.5 cm segment of muscle diaphragm containing the central tendon and the ribs was used.

1

Step 4.

rib was gently pinned to a holder fixed in the bottom of an organ bath filled with 5 ml of Tyrode's solution, while the central diaphragm tendon portion was attached to a PowerLab force transducer, which was maintained at 30°C and continuously gassed with 95% O₂/5% CO₂.

1

Step 5.

Isometric twitch contraction was elicited by electrical stimulation of muscle strips through silver electrodes, with 0.1 Hz frequency, 2 ms duration, and supramaximal voltage;

1

Step 6.

Tissue was rinsed four times with Tyrode's solution and incubated with 1 µM d-tubocurarine

1

Step 7.

After a 20 to 30-min stabilization, muscle length was readjust to give an optimal twitch tension, and 30 min later the effect of drugs was investigate.

1

Step 8.

Tyrode solution

Step 9.

mix in distilled water

☐ CONCENTRATION

11 Mass Percent Additional info: D-glucose

☐ CONCENTRATION

135 Mass Percent Additional info: NaCl

☐ CONCENTRATION

5 Mass Percent Additional info: KCl

☐ CONCENTRATION

0 Mass Percent Additional info: MgCl₆.H₂O

☐ CONCENTRATION

2 Mass Percent Additional info: CaCl₂

☐ CONCENTRATION

0 Mass Percent Additional info: Na₂HPO₄. H₂O

☐ CONCENTRATION

1 Mass Percent Additional info: NaHCO₃

Tyrode solution

Step 10.