

# Influenza A H3 virus TaqMan assay

Ian Mackay, Judy Northill

## Abstract

This assay is a modification to the World Health Organization's influenza A H3 TaqMan documented in 'WHO information for molecular diagnosis of influenza virus - update 1' (see file below). New primers were added and the WHO primers discarded.

**Citation:** Ian Mackay, Judy Northill Influenza A H3 virus TaqMan assay. **protocols.io**

dx.doi.org/10.17504/protocols.io.kmucu6w

**Published:** 08 Mar 2018

## Before start

- If using a different brand or model of real-time thermocycler, check that the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

## Protocol

### Oligonucleotides

#### Step 1.

Name	5'-3' SEQUENCE
H3hFor1	GGTACGGYTTTCAGGCAT
H3hRev1	TCAATCTGATGGAATTTCTCGTTG
H3h-1144dProbe	6FAM-CTGCTGCTTGTCCTCTCCCT-BHQ1

The oligoprobe is from the World Health Organization protocol linked below.

New primers were designed to improve assay performance.

🔗 LINK:

[http://www.who.int/entity/influenza/gisrs\\_laboratory/molecular\\_diagnosis\\_influenza\\_virus\\_humans\\_update\\_201403rev201505.pdf?ua=1](http://www.who.int/entity/influenza/gisrs_laboratory/molecular_diagnosis_influenza_virus_humans_update_201403rev201505.pdf?ua=1)

#### Step 2.

🧴 REAGENTS

SuperScript™ III Platinum™ One-Step qRT-PCR Kit [11732088](#) by [Life Technologies](#)

## Reaction Set-up

### Step 3.

The assay has been used on both a Rotor-Gene 6000 and a Rotor-Gene Q real-time thermocycler

Prepare sufficient mix for the number of reactions.

Include a suitable 'dead volume' as necessary if using a robotic dispenser.

#### MIX PREPARATION

Reagent	Vol (μL) X1	Final reaction concentration
Nuclease-free water	4.43	
H3hFor1 (200pmol/μL)	0.05	500nM
H3hRev1 (200pmol/μL)	0.05	500nM
H3h-1144dProbe (100pmol/μL)	0.03	150nM
2X Reaction Mix <sup>1</sup>	10.0	1X
ROX reference Dye (25μM) <sup>1,2</sup>	0.04	50nM
SuperScript <sup>TM</sup> III/Platinum <sup>TM</sup> Taq Mix <sup>1</sup>	0.4	
<b>TOTAL</b>	<b>15</b>	

<sup>1</sup>Superscript<sup>TM</sup> III Platinum<sup>TM</sup> One-step qRT-PCR kit; <sup>2</sup>See Guidelines

Dispense 15μL to each reaction well.

Add 5μL of template (extracted RNA, controls or NTC [nuclease-free water] ).

Total reaction volume is 20μL

## Amplification

### Step 4.

RT-PCR


50°C	5min	
95°C	2min	
95°C	3s	40X
60°C	30s*	

\*Florescence acquisition step

**Step 5.**

The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:

1. A **sigmoidal curve** – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
2. A **suitable level of fluorescence** intensity as measured in comparison to a positive control (y-axis)
3. A defined threshold ( $C_T$ ) value which the fluorescent curve has clearly exceeded (Fig.1 arrow) and which sits early in the log-linear phase and is  $<40$  cycles
4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a  $C_T$  value  $>40$  cycles is considered a negative result
5. NTCs should not produce a curve



**Figure 1.** Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle ( $C_T$ ) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.