

Transformation of E coli. with Heat Shock

Joshua Timmons

Abstract

Transformation of heat-shock competent E. coli cells. Adapted to fit the protocol followed by Northeastern_Boston

Citation: Joshua Timmons Transformation of E coli. with Heat Shock. **protocols.io**

dx.doi.org/10.17504/protocols.io.dqe5td

Published: 05 Sep 2015

Protocol

Step 1.

Thaw the appropriate amount of competent cells on ice.

Step 2.

Pre-chill the required number of empty 1.5 ml microcentrifuge tubes.

Step 3.

Pipet 50 µl aliquots of cells into the pre-chilled tubes.

Step 4.

Add 5-10 µl of a ligation reaction mix or 5 ng of pure plasmid DNA to each tube. Mix gently!

Step 5.

Incubate the tubes on ice for 30 min

⌚ DURATION

00:30:00

Step 6.

Heat shock the cells for 45 sec at 42°C

⌚ DURATION

00:00:45

Step 7.

Place the tubes immediately on ice for at least 2 min

⌚ DURATION

00:02:00

Step 8.

Add 800 µl of SOC medium to each tube and incubate for 1 hour at 37°C

⌚ DURATION

01:00:00

Step 9.

Transfer the cultures to 1.5 ml microcentrifuge tubes and spin for 1 min at 6000 rpm.

⌚ DURATION

00:01:00

Step 10.

Remove 800 µl of the supernatant and resuspend the pellet.

Step 11.

Plate out the suspension on a LB agar plate containing the appropriate antibiotic.

Step 12.

Incubate the plates overnight at 37°C.