

Capped RNA E2050 Mixture

New England Biolabs

Abstract

Citation: New England Biolabs Capped RNA E2050 Mixture. **protocols.io**

dx.doi.org/10.17504/protocols.io.crzv75

Published: 04 Feb 2015

Protocol

Step 1.

Nuclease-free water so that total volume is **20 µl**

Step 2.

NTP Buffer Mix, **2 µl** (2 mM each NTP final)

 **AMOUNT**

2 µl Additional info:

Step 3.

Cap Analog (40 mM), **4 µl** (8 mM final)

 **AMOUNT**

4 µl Additional info:

Step 4.

Template DNA, X µl (1 µg)

Step 5.

T7 RNA Polymerase Mix, **2 µl**

 **AMOUNT**

2 µl Additional info: