

# Adult mouse pancreas cell dissociation (on ice)

#### **Andrew Potter**

#### **Abstract**

This procedure is used to dissociate adult (10 wk.) mouse pancreas into single cells. The procedure is carried out on ice in order to maintain a more authentic gene expression profile. It is a two-layered dissociation, with each layer consisting of 5 mg/mL type 4 collagenase. The yield is 4400 cells/mg with 94% viability.

Citation: Andrew Potter Adult mouse pancreas cell dissociation (on ice). protocols.io

dx.doi.org/10.17504/protocols.io.r8rd9v6

Published: 31 Jul 2018

# **Guidelines**

# Type 4 Collagenase Enzyme Mix (2 x 1 mL)

50 μL Type 4 Collagenase 100 mg/mL (5 mg/mL final conc.)

5 μL DNAse (125 U)

5 μL 1 M CaCl<sub>2</sub> (5 mM final conc.)

10 μL 10% BSA/PBS (0.1% BSA final conc.)

930 µL DPBS

# Reagents

Collagenase Type 4 - Worthington (LS004186) - 100 µL aliquots of 100 mg/mL, frozen at -80 °C

DNAse 1 - Applichem (A3778) - 10 µL aliquots each with 250 U, frozen at -80 °C

DPBS - ThermoFisher (cat. #14190)

Red Blood Cell Lysis Buffer - Sigma (R7757) Trypan Blue Solution 0.4% - Gibco (15250061)

# **Required Equipment & Consumables:**

Refrigerated centrifuge

Pipettes and pipet tips (MLS) 15, 50 ml Conicals (MLS) 1.5 mL tubes (MLS) 30 µM filters (MACS SmartStrainers, 130-098-458)

Petri dishes (MLS)
Razor blades (MLS)
Ice bucket w/ice (MLS)
Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

## The protocol workflow is as follows:

- 1. Isolate pancreas
- 2. First layer
- 3. Second layer
- 4. Preparing cells for Chromium/DropSeq

### **Before start**

- -Prepare enzyme mixes and leave on ice.
- -Cool centrifuges to 4 °C.
- -Isolate and transport tissue in ice-cold DPBS.

# **Protocol**

## Step 1.

Dissect pancreas and place in ice-cold PBS.

#### Step 2.

Mince tissue thoroughly on petri dish on ice (2 min) until fine paste.

**O** DURATION

00:02:00: mince on ice

# Step 3.

Weigh out 18 mg tissue and add to tube with 1 mL Type 4 collagenase enzyme mix.

**AMOUNT** 

18 mg: minced pancreas tissue

## Step 4.

Incubate on ice. Shake vigorously every 30 seconds for the first two min to re-suspend tissue.

O DURATION

00:00:30 : shake vigorously for first 2 min

#### Step 5.

At two mins, begin triturating 10x every min. Continue triturating on ice for 20 min.

**O DURATION** 

00:20:00: triturate on ice

#### Step 6.

After incubating 20 min, let chunks settle for 1 min on ice.

**O** DURATION

00:01:00: let chunks settle

#### Step 7.

Pipet top 75% (750  $\mu$ L) of supernatant containing released cells onto a 30  $\mu$ M filter on a 50 mL conical, on ice.

**■** AMOUNT

750 µl : save supernatant

# Step 8.

Rinse filter with 5 mL ice-cold PBS/BSA 0.04%. Save filter and flow-through for next steps.

**AMOUNT** 

5 ml: ice-cold PBS/BSA 0.04%

#### Step 9.

To residual tissue chunks add additional 1 mL type 4 collagenase enzyme mix.

**■** AMOUNT

1 ml: type 4 collagenase mix

#### Step 10.

Continue triturating on ice 10x every min for 30 additional min (50 min total digest time).

**O DURATION** 

00:30:00: triturate on ice

#### Step 11.

Triturate and add entire volume to same 30  $\mu$ M filter on 50 mL conical. Rinse filter w/5 mL ice-cold PBS/BSA 0.04%.

**■** AMOUNT

5 ml: ice-cold PBS/BSA 0.04%

# Step 12.

Transfer flow-through to 15 mL conical. Spin 300 g for 5 min at 4 °C.

**O DURATION** 

00:05:00 : spin 300 g at 4 °C

#### Step 13.

Remove supernatant and re-suspend in 100 µL ice-cold PBS/BSA 0.04%.

#### Step 14.

Add 900 µL RBC lysis buffer. Triturate 20x and incubate 2 min. on ice.

**O DURATION** 

00:02:00 : incubate on ice

Step 15.

Add 10 mL ice-cold PBS/BSA 0.04% to dilute RBC lysis buffer.

Step 16.

Spin 300 g for 5 min at 4 °C. Remove supernatant.

**O DURATION** 

00:05:00 : spin 300 g

Step 17.

Re-suspend in 100  $\mu$ L ice-cold PBS/BSA 0.04%. Analyze viability and yield using a hemocytometer with trypan blue. Adjust concentration to 1,000 cells/ $\mu$ L for Chromium or 1,000 cells/ $\mu$ L for DropSeq.

**■** AMOUNT

 $100~\mu l$  : ice-cold PBS/BSA 0.04%