



## 16S Illumina Primer Plate Resuspension (100 µM)

Rodolfo Salido Benitez<sup>1</sup>

<sup>1</sup>Knight Lab @ UC San Diego Health

[dx.doi.org/10.17504/protocols.io.ukgeutw](https://doi.org/10.17504/protocols.io.ukgeutw)

 Rodolfo Salido Benitez  
Knight Lab @ UC San Diego Heal... 

### ABSTRACT

The following protocol describes resuspension of lyophilized primers from (4) 96-Well plates shipped from Integrated DNA Technologies (IDT) at 3nmol per well using the epMotion 5075.

For information about ordering consult the [Earth Microbiome Projects website](#).

Barcoded primer constructs for 16S V4 region are also described in the [Earth Microbiome Project website](#).

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
ep T.I.P.S. Motion Racks 20 - 300 µL w/ filter	0030014456	Eppendorf
PCR Clean Water	<a href="#">View</a>	
IDT Lyophilized Primer Plate	<a href="#">View</a>	Integrated DNA Technologies
epMotion Reservoir 30 mL	960051009	Eppendorf

### MATERIALS TEXT

- (4) ep T.I.P.S. Motion Racks 20 -300 µL w/ filter
- (4) IDT 96-Well Lyophilized Primer Plates (3 nmoles of primer per well)
- (1) epMotion Reservoir 30 mL

### BEFORE STARTING

Please wear at least the minimum required personal protective equipment.

Ensure that all necessary kit components are available as well as user-supplied consumables.

Remove nuclease and nucleotide contamination from work surfaces and instruments prior to starting using an appropriate solution, such as RNase AWAY™ (Thermo Scientific™ catalogue: 700511), followed by wiping with 70% to 100% molecular biology grade ethanol to remove additional contaminants.

#### Prepare reagents

- 1 In a sterile 30mL reservoir, add 11.725 ml (required minimum volume) of PCR Clean Water to resuspend (4) 96-Well plates.

 11.725 ml PCR Clean Water

#### Prepare Plates

- 2 Centrifuge all 96-Well primer plates to ensure lyophilized sample is at the bottom of the well. (3nmol primer per well)


3

**EQUIPMENT**

epMotion 5075  
Liquid Handling  
Eppendorf 5075000962

Follow the diagram below while setting up the epMotion worktable.

Worktable



Place (4) boxes of ep.T.I.P.S. Motion Racks 20 - 300  $\mu$ L w/ filter in slots A2-A5

Place (4) 96-Well plates of lyophilized primers from IDT in slots B2-B5

Place 30 mL reservoir with PCR Clean Water in slot 1 of the Reservoir Rack and place Rack in slot C5

## Execute automated protocol

4 Remove box lids and plate foils and execute protocol.

(Protocol must be imported to epBlue software prior to attempting to execute it. *epBlue 40.6 or later*)

☐ **Application\_4\_Resus\_100uM\_Primers\_181003\_105824.export6**

The automated protocol transfers 30 $\mu$ L of PCR Clean Water into the (4) different 96-Well Primer Plates and mix pipettes for 5 cycles to resuspend lyophilized primers to a 100  $\mu$ M

## Seal and store plates.

5 Remove plates from worktable and seal with storage aluminum foils.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited