



Nov 01,  
2019

# Enzymatic Spectrophotometry Assay to Measure D-, L-, and Total Lactate in Canine Feces [↗](#)

PLOS One

Amanda Blake<sup>1</sup>, Jonathan Lidbury<sup>1</sup>, Joerg Steiner<sup>1</sup>, Jan Suchodolski<sup>1</sup>

<sup>1</sup>Texas A&M University - College Station

**1** Works for me [dx.doi.org/10.17504/protocols.io.hrbq55w](https://doi.org/10.17504/protocols.io.hrbq55w)



Amanda Blake  
Texas A&M University - College Station



## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0224454>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Blake AB, Guard BC, Honneffer JB, Lidbury JA, Steiner JM, Suchodolski JS (2019) Altered microbiota, fecal lactate, and fecal bile acids in dogs with gastrointestinal disease. PLoS ONE 14(10): e0224454. doi: [10.1371/journal.pone.0224454](https://doi.org/10.1371/journal.pone.0224454)

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
Triethanolamine	T58300	Sigma Aldrich
Trichloroacetic acid solution	T0699	Sigma Aldrich
D-/L-Lactic acid kit	11112821035	

### Aliquot feces

- 1 Aliquot 120-130 mg feces into 1.5 mL microcentrifuge tube. Store at -80°C until deproteinization.

### Begin deproteinization. Add TEA buffer

- 2 Make 0.1 M triethanolamine buffer (pH 9.15) using scale for accuracy (to the 0.01 mg). Then add 750 µL buffer to each aliquot. Vortex and place in refrigerator (-4°C) for three hours, vortexing every hour to ensure thorough mixing.



Triethanolamine

by Sigma Aldrich

Catalog #: T58300

### Centrifuge #1

- 3 Centrifuge samples at 13,000 x g for 5 minutes at 4°C. Carefully pipette 495 µL of supernatant into a new 1.5 mL microcentrifuge tube. Store supernates at -80°C overnight to provide ample time for next processing steps.



From empirical experience, the extra freeze-thaw cycles favorably increased protein pellet size while not affecting lactate concentration.

#### Add TCA

- 4 Thaw samples at room temperature for 20-40 minutes. Once thawed, add 10  $\mu$ L of 6 M trichloroacetic acid to each sample. Vortex for 10 seconds and place immediately in ice bath for 20 minutes.



Performed this step and following steps (up to dilution with TEA) under a chemical fume hood with proper personal protective equipment.

<http://www.sigmaaldrich.com/MSDS/MSDS/DisplayMSDSPage.do?country=US&language=en&productNumber=T0699&brand=SIGMA&PageToGoToURL=http%3A%2F%2Fwww.sigmaaldrich.com%2Fcatalog%2Fproduct%2Fsigma%2Ft0699%3Flang%3Den>



#### Trichloroacetic acid solution

by Sigma Aldrich

Catalog #: T0699

#### Centrifuge #2

- 5 Vortex samples again upon removing from ice bath, and centrifuge at 4,500 x g for 20 minutes at 4°C. After centrifugation, a protein pellet was noted in the bottom of each tube. Next, remove 400  $\mu$ L of supernatant and place in a clean 2 mL microcentrifuge tube.

#### Dilute samples

- 6 Add 1600  $\mu$ L 0.1 M triethanolamine buffer (pH 9.15) to each sample to achieve a neutral or alkaline pH (between 7 and 10). These deproteinized fecal extracts were either used immediately for lactate analysis or stored in -80°C for later use.

#### Make standards

- 7 Make standard dilutions with ultra-pure water and D- or L-lactate provided in enzymatic kit (D-/L-Lactate Enzymatic Kit, R-Biopharm Inc.)

Standard	Dilution	Concentration (g/L)
1	1:4	0.05175
2	1:10	0.0207
3	1:20	0.01035
4	1:40	0.005175
5	1:80	0.0025875
Stock concentration = 0.207 g/L		



#### D-/L-Lactic acid kit

Catalog #: 11112821035

### Make Master Mix

- 8 Make master mix (MM) solution on an as needed basis depending on the number of samples. In a 15 or 50 mL tube, mix together (100 x n)  $\mu$ L solution 1, (20 x n)  $\mu$ L solution 2, and (2 x n)  $\mu$ L solution 3, with n equal to the total number of blanks, standards, and unknown samples.

Solution Number	Description*
1	approx. 30 ml solution, consisting of: glycylglycine buffer, pH approx. 10.0; L-glutamic acid, approx. 440 mg
2	approx. 210 mg NAD, lyophilizate, reconstituted in 6 ml redist. Water
3	approx. 0.7 ml glutamate-pyruvate transaminase suspension, approx. 1100 U
4	approx. 0.7 ml D-lactate dehydrogenase solution, approx. 3800 U
5	approx. 0.7 ml L-lactate dehydrogenase solution, approx. 3800 U
*Descriptions obtained from manufacturer protocol	

### Plate your samples and incubate

- 9 Pipette blanks, standards, and samples in duplicate onto 96-well plate and place on plate shaker for one minute. Shake at low enough speed that a plate cover is not needed. Then incubate the plate in plate incubator at 25°C for 15 minutes.

Well type	Volumes
blanks	122 $\mu$ L MM, 100 $\mu$ L water
standards (D- and L-lactate)	122 $\mu$ L MM, 100 $\mu$ L standards 1-5
unknown samples	122 $\mu$ L MM, 100 $\mu$ L fecal extract
MM = master mix	

### Plate reader and settings

- 10 Use path length adjustment settings in plate reader software to account for using a protocol originally designed for use with 1 cm diameter cuvettes. Set standards curve to quadratic and adjust lactate concentrations by subtracting the absorbance difference of the blank.



#### Gen5 2.07

Windows

by BioTek® Instruments, Inc.

### A1 and D-LDH

- 11 After incubation, take the first absorbance reading (A1) at 340 nm wavelength. Then immediately add 2  $\mu$ L D-LDH to each well, shake plate for one minute, and then incubate at 25°C for 30 minutes.

### A2 and L-LDH

- 12 Take second absorbance reading (A2) at 340 nm wavelength. Add 2  $\mu$ L L-LDH to each well, shake for one minute, and incubate at 25°C for 30 minutes.

- 13 Take third absorbance reading at 340 nm wavelength.

#### Data processing

- 14 Calculate standard curve and lactate concentrations using software above. Then, if D- or L-lactate concentrations are below their respective lower limits of quantification, adjust concentrations to their respective cutoff values. Ours were 0.002 g/L for D-lactate and 0.0007 g/L for L-lactate. Also adjust final lactate concentrations based on starting weight of feces, dilution factor, and dry matter content.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited