

Immunoprecipitation using Protein A/G Magnetic Beads

New England Biolabs

Abstract

This is the protocol for immunoprecipitation using Protein A/G Magnetic Beads.

Citation: New England Biolabs Immunoprecipitation using Protein A/G Magnetic Beads. protocols.io

dx.doi.org/10.17504/protocols.io.crav2d

Published: 02 Feb 2015

Before start

Use 25 μ l of Protein A/G Magnetic Beads per 200 μ l of crude cell lysate containing 200-500 μ g of total protein in a standard immunoprecipitation protocol. It is important to increase the volume of beads proportionately for larger cell lysate volumes.

Materials

- Frotein A Magnetic Beads 1 ml S1425S by New England Biolabs
- Frotein G Magnetic Beads 1 ml S1430S by New England Biolabs

Protocol

Cell Lysis

Step 1.

Rinse a 60 mm culture dish of confluent cells with PBS

Cell Lysis

Step 2.

Lyse the cells with 0.5 ml cold Immunoprecipitation Buffer

AMOUNT

1 ml Additional info:

₹ PROTOCOL

. NEB Immunoprecip. Buffer

CONTACT: New England Biolabs

Step 2.1. 50 mM NaCl

Step 2.2.

10 mM Tris-HCl (pH 7.4)

Step 2.3.

1 mM EDTA

Step 2.4.

1 mM EGTA (pH 8.0)

Step 2.5.

0.2 mM sodium ortho-vanadate

Step 2.6.

0.2 mM PMSF

Step 2.7.

1% Triton X-100

Step 2.8.

0.5% NP-40

Cell Lysis

Step 3.

Maintain constant agitation for 30 minutes at 4°C

O DURATION

00:30:00

Cell Lysis

Step 4.

Scrape the cells from the dish

Cell Lysis

Step 5.

Sonicate on ice for 5 seconds; repeat 4 times

Cell Lysis

Step 6.

Centrifuge for 5 minutes at 4°C

© DURATION

00:05:00

Cell Lysis

Step 7.

Assay for total protein then adjust concentration to approximately 1 mg/ml with Immunoprecipitation Buffer

NOTES

New England Biolabs 31 Jan 2015

The supernatant is the crude cell lysate.

Immunoprecipitation

Step 8.

In a 1.5 ml microcentrifuge tube, add $25~\mu l$ Protein A/G Magnetic Beads to $200~\mu l$ of crude cell extract

NOTES

New England Biolabs 01 Feb 2015

Steps 8-12 pre-clear crude cell extract of proteins which can bind non-specifically to the beads

Immunoprecipitation

Step 9.

Gently vortex

Immunoprecipitation

Step 10.

Incubate at 4°C for 1 hour

© DURATION

01:00:00

Immunoprecipitation

Step 11.

Apply magnetic field for 30 seconds to pull beads to the side of the tube

O DURATION

00:00:30

Immunoprecipitation

Step 12.

Pipette supernatant to a clean 1.5 ml microcentrifuge tube and discard the beads.

Immunoprecipitation

Step 13.

Add 1-5 μg of antibody to crude cell lysate

Immunoprecipitation

Step 14.

Vortex

Immunoprecipitation

Step 15.

Incubate at 4°C for 1 hour

© DURATION

01:00:00

P NOTES

New England Biolabs 01 Feb 2015

If monoclonal antibodies are used, add 5 μ g rabbit anti-mouse IgG antibody. Vortex and incubate an additional 30 minutes at 4°C. Alternatively, Protein G Magnetic Beads (NEB #S1430S) can be used for immunoprecipitations with monoclonal antibodies.

Immunoprecipitation

Step 16.

Add 25 µl of Protein A/G Magnetic Beads suspension

Immunoprecipitation

Step 17.

Gently vortex

Immunoprecipitation

Step 18.

Incubate with agitation for 1 hour at 4°C.

© DURATION

01:00:00

Immunoprecipitation

Step 19.

Apply magnetic field to pull beads to the side of the tube

Immunoprecipitation

Step 20.

Carefully pipette to remove supernatant

Immunoprecipitation

Step 21.

(wash #1) Wash with **500 µl** of Immunoprecipitation Buffer by gentle vortex

Immunoprecipitation

Step 22.

(wash #1) Apply magnetic field then remove supernatant and discard

Immunoprecipitation

Step 23.

(wash #2) Wash with **500 µl** of Immunoprecipitation Buffer by gentle vortex

Immunoprecipitation

Step 24.

(wash #2) Apply magnetic field then remove supernatant and discard

Immunoprecipitation

Step 25.

(wash #3) Wash with 500 µl of Immunoprecipitation Buffer by gentle vortex

Immunoprecipitation

Step 26.

(wash #3) Apply magnetic field then remove supernatant and discard

Immunoprecipitation

Step 27.

Resuspend bead pellet in **30 µl** of 3X SDS Sample Loading Buffer

■ AMOUNT

30 µl Additional info:

O DURATION

00:05:00

PROTOCOL

. NEB 3X SDS Sample Loading Buffer

CONTACT: New England Biolabs

Step 27.1.

187.5 mM Tris-HCl (pH 6.8)

Step 27.2.

6% (w/v) SDS

Step 27.3.

30% glycerol

Step 27.4.

150 mM DTT

Step 27.5.

0.03% (w/v) bromophenol blue

Step 27.6.

2% β-mercaptoethanol

Immunoprecipitation

Step 28.

Incubate sample at 70°C for 5 minutes.

O DURATION

00:05:00

Immunoprecipitation

Step 29.

Apply magnetic field to sample then load supernatant on SDS-PAGE gel and electrophorese.