



Nov 27,  
2019

## Western Blot

Athan Dial<sup>1</sup>

<sup>1</sup>McMaster University

1

Works for me

[dx.doi.org/10.17504/protocols.io.rh3d38n](https://doi.org/10.17504/protocols.io.rh3d38n)



Enoch Tin ⚡

### ABSTRACT

Biochemical technique to separate protein based on molecular weight. Assumed to have homogenized samples.

### SDS-PAGE

- 1 Pull out comb gently and straight, avoid collapsing wells
- 2 Follow instructions provided with power pack to set up gels
  - Short plate facing in on apparatus
  - Press against gasket
  - Press down and close arms simultaneously
  - Put unit into bucket
  - If not using 2 gels, use thick plast plate as the "dummy plate"
- 3 Fill inner chamber with 'fresh' 1x running buffer high enough to cover smaller plate and enter wells
- 4 Load protein ladder in right most lane  
\*Note: gels are reversed, final product will have ladder on leftmost lane
- 5 Load samples with gel-loading tips at consistent amount  
  
Tips: avoid bubbles and draw/eject samples slowly
- 6 Pour "used" 1x running buffer to fill bucket
- 7 Place the lid on the tank with matched colour-coded leads. Run the machine until protein ladder sufficiently separated.

⚡ 25 °C Room temperature

🕒 02:15:00 Constant 60 V

## Prep Transfer Buffer (1L per container)

- 8 Make 1 L solution in "1x transfer buffer with MeOH" bottle

🧴 700 ml Distilled Water

🧴 200 ml Methanol

🧴 100 ml 10x Transfer Buffer

- 9 Put transfer buffer into -20°C freezer #2 (bottom rack)

## Transfer

- 10 Fill wide plastic tray with 1x transfer buffer

- 11 Soak filter papers (6 in total per cassette) in 1x transfer buffer in tray

- 12 Carefully remove gel from between glass plates. Use the green plastic wedge.

- 13 Using tweezers, obtain PVDF membrane and place in methanol on rocker. This activates the membrane.

- 14 Soak gels in transfer buffer for 5 minutes to equilibrate. Soak sponges as well.

- 15 Assemble immunoblot 'sandwich' in following order:

- Blot cassette with black (cathode) side down
- 2 sponges
- 3 sheets of soaked filter paper (roll out bubbles)
- Gel (with ladder on the right, dye front on the bottom)
- With tweezers, place PVDF membrane over gel
- 3 sheets of soaked filter paper (roll out bubbles)
- 2 sponges

Close the cassette. Leave the cassette in transfer buffer until all cassettes are assembled

- 16 Set up the power pack.

- Fill the bucket with the transfer buffer (should be about 3/4 full)
- Put in magnetic stir bar
- Put in ice pack into the bucket
- Ensure electrodes are in correct placement

Start reaction in the cold room. Small bubbles should be seen when the reaction starts

🕒 01:30:00 Constant 90V

## Block

- 17 Block reagents and concentration will vary depending on antibodies.  
Usually 10 mL of block solution per plate.

- 18 Place membrane onto plate with tweezers. Add used amido black. Place plate on the rocker for 1 minute.

19 Draw up excess amido black back into the bottle. Rinse the membrane with distilled water.

20 Cut the membrane with exacto-knife (if necessary). Take a picture.

21 Place membrane on rocker with TBS-T to keep the membrane moist.

22 Remove previous TBS-T. Add block solution and incubate on rocker

🕒 01:00:00 Incubate on rocker

#### Primary Antibody Incubation

23 Defrost primary antibody solution 5 mins before completion of block incubation

24 Remove block solution and add primary antibody solution.

25 Incubate on a rocker within the cold room **overnight**. Keep the bottle with the membrane for reference purposes.

#### Next Day

26 ---

🕒 14:00:00 Incubate with primary antibody (min)

🕒 18:00:00 Incubate with primary antibody (max)

#### Secondary Antibody Incubation

27 Remove primary antibody with 10mL pipette. Keep the solution and put it back into its respective tube.

28 Wash membrane:

Add 1x TBS-T to the membrane and incubate on the rocker

🕒 00:05:00 Wash/incubate on rocker


29 Perform wash step 3 times in total

30 Prepare secondary antibody solution. Diluent and % mass will vary depending on primary antibody used.  
10 mL solution per plate.

Do not contaminate antibodies.




31 Add secondary antibody to the membrane

🕒 01:00:00 Incubate on rocker


- 32 Wash membrane:  
Add 1x TBS-T to the membrane and incubate on the rocker  
 **00:05:00 Wash/Incubate on rocker**

- 33 Perform wash step 3 times in total

#### ECL Detection

- 34 Keep membranes in TBS-T but remove from the rocker while preparing the ECL solution.
- 35 Set-up the imager.
  - Turn on the computer
  - Autofocus with the Mac card
- 36 Clean the square plate with EtOH. Ensure no particular matter
- 37 Prepare ECL solution in a black tube. Mix using the pipette.  
**Use different pipette tips to prevent contamination.**  
 **750 µl Luminol**  
 **750 µl Clarity**
- 38 Dap off excess TBS-T with KimWipe. Place membrane in square plate. Pipette out some of the ECL solution. In drop-wise fashion, just cover the membrane in the ECL solution.
- 39 Place container over the plate to avoid light  
 **00:05:00 Incubate in solution**
- 40 Using tweezers, put the membrane on the glass. Place the glass on the 3rd level in the machine in the middle of the tray.
- 41 Take preview image to see the ladder. Save the image.
- 42 Auto-expose and save the image

#### Membrane Storage

- 43 Soak membrane in 100% methanol and leave on a KimWipe to dry. Wrap the membrane in plastic wrap and label it.  
 **25 °C Store at RT**



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited