

Lethal Phase - Immunohistochemistry

Sonia M. Hall

Abstract

Experiment was conducted to determine the stage at which mutant embryos arrest their development. All *Drosophila melanogaster* lines were crossed to balancers with a YFP marker for unambiguous identification of mutant and non-mutant embryos. Embryos were collected on apple juice agar plates with a smear of yeast paste for 1hr and aged for 17 hours at 25°C. Embryos were fixed and stained as described by Fehon et al., 1991 using a modified rotational speed of 150rpm.

Citation: Sonia M. Hall Lethal Phase - Immunohistochemistry. **protocols.io**

dx.doi.org/10.17504/protocols.io.cxfxjm

Published: 29 Jul 2015

Protocol

Step 1.

Wash embryos and place into wash basket. Dechorionate by pouring bleach directly onto embryos. Let sit for three minutes.

 DURATION

00:03:00

Step 2.

Rinse embryos in basket with ddH₂O until bleach is completely washed away.

Step 3.

Place embryos into glass vial containing 6ml of Heptane and 6ml of 0.4% Paraformaldehyde. Shake on orbital rotator at 150rpm for 20 minutes.

 DURATION

00:20:00

Step 4.

Upon completion of 20 minute fixation, allow heptane and paraformaldehyde to separate. Embryos will be trapped at the interface between the two solutions. Use a glass pipette to remove the bottom (paraformaldehyde) layer be sure to not disrupt the embryos at the interface.