

# untitled protocol

Lutz Berwanger

## Abstract

In vitro transcriptions assay for cyanobacterial RNA polymerase. This protocol can be used either to do an in vitro transcription with small DNA fragments with desired promoter region or to use it with a promoter region integrated in a plasmid.

**Citation:** Lutz Berwanger untitled protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.mpnc5me

**Published:** 17 Jan 2018

## Guidelines

RNA experiments! Always wear gloves and prepare probes under RNase free and cooled conditions!

## Before start

Prepare all Buffer, measure RNAP and DNA concentrations

## Materials

 in vitro transcription buffer by Contributed by users

## Protocol

### Preparations

#### Step 1.

*In vitro* transcription Buffer:

20 mM Tris-HCl pH 7,9

40 mM KCl,

15 mM MgCl<sub>2</sub>

## Prepartions

### Step 2.

Polyacrylamid Gel:

10-20% 19:1 acrylamide:bisacrylamide

8M urea

1x TBE

### NOTES

**Lutz Berwanger** 17 Jan 2018

For smaller fragments use higher concentrations of 19:1 acrylamide:bisacrylamide.

## Prepartions

### Step 3.

Loading dye should contain at least 50% formamide or 8M urea. You can use commercial RNA loading dye, also.

Example recipie:

8M urea saturated with formamide

10mM EDTA

1xTBE,

bromophenol blue and xylene cyanol (amounts are empirical, just to make it visible on the gel)

## In vitro Transcription

### Step 4.

Mix:

100 nM RNAP

10 nM promoter-containing DNA

100µM NTPs

#### 📌 NOTES

**Lutz Berwanger** 17 Jan 2018

If desired one could use P32 labelled NTPs. Then it is recommendable to decrease the concentration of the same cold NTP to 10uM, to increase labelling efficiency.

### In vitro Transcription

#### Step 5.

Incubate the reaction at 34°C for at least 3 hours.

Stop with equal volume of the loading dye. Before loading onto denaturing PAA gel boil the Samples for 5 min at 95°C and let them cool on ice.

#### 📌 NOTES

**Lutz Berwanger** 17 Jan 2018

If desired, one could expand reaction duration.

### Visualization

#### Step 6.

Visualization of the products on the gel, can be done with nucleic acid gel stain of your choice. For small products it is recommended to stain the gel over night.

#### 🧴 REAGENTS

GelRed™ Nucleic Acid Gel Stain, 10,000X in Water [G-725](#) by [Gold Biotechnology](#)

✓ SYBR Gold Nucleic Acid Gel Stain [S-11494](#) by Contributed by users

#### 📌 NOTES

**Lutz Berwanger** 17 Jan 2018

Two common used stains are indicated.

## Warnings

Harmless