

# An. gambiae complex species identification by PCR

Fabio Gomes

## Abstract

Scott, JA, WG brogdon, and FH Collins, *Identification of single specimens of the Anopholes gambiae complex by the polymerase chain reaction*. Am J Trop Med Hyg, 1993. **49**(4): p. 520-9

**Citation:** Fabio Gomes An. gambiae complex species identification by PCR. **protocols.io**

dx.doi.org/10.17504/protocols.io.d7i9kd

**Published:** 17 Nov 2015

## Guidelines

Guideline for Identification

M Form: 370bp

S Form: 260bp

## Before start

### PRIMERS - 5' TO 3'

									Size of fragment (BP)	Tm
Universal (forward)	GTG	TGC	CCC	TTC	CTC	GAT	GT			<b>58.3</b>
GA - gambiae (reverse)	CTG	GTT	TGG	TCG	GCA	CGT	TT		<b>390</b>	<b>59.3</b>
ME - merus & melus (reverse)	TGA	CCA	ACC	CAC	TCC	CTT	GA		<b>466/464</b>	<b>57.2</b>
AR - arabiensis (reverse)	AAG	TGT	CCT	TCT	CCA	TCC	TA		<b>315</b>	<b>47.4</b>
QD - quadriannulatus (reverse)	CAG	ACC	AAG	ATG	GTT	AGT	AT		<b>153</b>	<b>42.7</b>

## Materials

🐛 HhaI - 2,000 units [R0139S](#) by [New England Biolabs](#)

Platinum PCR SuperMix 11306-016 by [Thermo Scientific](#)

## Protocol

### PCR

#### Step 1.

Assemble Mix

Reagent	Per rxn:
Platinum Supermix	16.5
UN (20uM)	0.5
GA (20 uM)	0.5
AR (20uM)	0.5
QD (20uM)	0.5
ME (20uM)	0.5
Template	1.0
<b>Total</b>	20ul

## PCR

### Step 2.

Run PCR setup

**Hot Start: 94° for 2 min**

**30 cycles of: 94° for 30 sec**

**50° for 30 sec**

**72° for 30 sec**

**Final extension: 72° for 5 min**

**4° forever**

## PCR

### Step 3.

Save 5uL for running a gel

## Restriction Digestion

### Step 4.

Take 10-15 uL PCR product, complete to 17.5 uL with NFW.

Add 2 uL 10x HhaI Buffer and 0.5 uL HhaI

Incubate at 37oC Overnihgt

## Restriction Digestion

### Step 5.

Run Agarose Gel using 100 bp ladder