

# Modified KAPA Express Extract

Michael Crone

## Abstract

KAPA Express Extract is a novel thermostable protease and buffer system that allows for the extraction of PCR-ready DNA from various tissue types in as little as 10 min. This protocol uses half the normal reaction mix.

**Citation:** Michael Crone Modified KAPA Express Extract. **protocols.io**

dx.doi.org/10.17504/protocols.io.czw7d

**Published:** 14 Jun 2015

## Guidelines

When running program on PCR machine incubate at 75°C and then 95°C consecutively and then cool the reaction products to 4°C.

## Protocol

### Reaction Setup

#### Step 1.

#### Reaction Setup

	Volume
KAPA Express Extract Buffer	5 µl
1U/µl KAPA Express Extract Enzyme	1 µl
PCR-Grade Water	44 µl
Sample	

### ✓ PROTOCOL

#### . [KAPA Express Extract Reaction Mix](#)

CONTACT: [Michael Crone](#)

#### Step 1.1.

10 X KAPA Express Extract Buffer, **5 µl**

📄 [AMOUNT](#)

5 µl Additional info:

#### Step 1.2.

1U/µl KAPA Express Extract Enzyme, **1 µl**

📄 [AMOUNT](#)

1 µl Additional info:

**Step 1.3.**

Water, **44 µl**

 **AMOUNT**

44 µl Additional info:

**Step 1.4.**

Sample

**Lysis****Step 2.**

Incubate in a thermocycler for 10 minutes at 75°C. During this step, cells are lysed, nucleases and proteins degraded and DNA released.

 **DURATION**

00:10:00

**Inactivation****Step 3.**

Incubate for 5 minutes at 95°C to inactivate the thermostable KAPA Express Extract protease.

 **DURATION**

00:05:00

**Centrifugation****Step 4.**

Vortex reaction product for 2-3 sec. Centrifuge at high speed for 1 minute to pellet debris.

Transfer DNA-containing supernatant to a fresh tube.

Use 1 µl of DNA extract directly in a 25 µl PCR, without quantification.

Dilute in TE Buffer for long-term storage at -20°C.