# Multiplex 1 PCR-SSP - CR1(rs3737002; rs11118131; rs11118167; rs17047660)

# Gabriela Canalli Kretzschmar, Luana Caroline Oliveira, Angelica Beate Winter Boldt

#### **Abstract**

This is a multiplex PCR-SSP for amplification of the follow SNPs of CR1 gene:

rs3737002;

rs11118131;

rs11118167;

rs17047660

Citation: Gabriela Canalli Kretzschmar, Luana Caroline Oliveira, Angelica Beate Winter Boldt Multiplex 1 PCR-SSP -

CR1(rs3737002; rs11118131; rs11118167; rs17047660). protocols.io

https://www.protocols.io/view/multiplex-1-pcr-ssp-cr1-rs3737002-rs11118131-p49dqz6

Published: 18 Jun 2018

#### **Before start**

- 1- Wear clean gloves;
- 2- Clean pipettes and stand with hypochlorite and 70% alcohol;
- 3- Defreeze DNA samples and reagents (exception: Tag polymerase);
- 4- Gently mix the DNA samples and pass the reagents briefly on the vortex.
- 5- Centrifuge the DNA samples and reagents with a spin;
- 6- Keep reagents on ice;
- 7- Leave microtubes ready for preparation of the mix, keep the same ones identified and on ice.
- 8- Leave application plates properly identified.

# **Materials**

- 1x Coral Buffer by Qiagen
- ✓ dNTP mix by Contributed by users
- custom made primers by Contributed by users
- ✓ 20 ng of genomic DNA by Contributed by users
- ✓ Ultrapure Water by Contributed by users
- Tag DNA Polymerase Platinum by <u>Invitrogen Thermo Fisher</u>

#### **Protocol**

#### Step 1.

The amounts described in this protocol are for one sample. Multiply them for the number of samples to be investigated.

#### Remember:

- Each sample shall be tested in all four mix reactions 1, 2, 3 and 4 (4 tubes).
- Always make plus 10% mix.

This means that each sample shall be distributed in 4 microtubes (for 4 reactions).

- Label properly the microtubes for the mixes, identifying them:
  - MIX 1 (in28-rs11118167 C/ ex29-rs17047660 A; ex26-rs3737002 C/ in26-rs11118131 C)
  - MIX 2 (in28-rs11118167 T/ ex29-rs17047660 A; ex26-rs3737002 T/ in26-rs11118131 T)
  - MIX 3 (in28-rs11118167 T/ ex29-rs17047660 G; ex26-rs3737002 C/ in26-rs11118131 T)
  - MIX 4 (in28-rs11118167 C/ ex29-rs17047660 G; ex26-rs3737002 T/ in26-rs11118131 C)
- Keep these microtubes in ice.

The sizes of the fragments are:

in28-rs11118167 / ex29-rs17047660 = 746 bpex26-rs3737002 / in26-rs11118131 = 457 bp

# Step 2.

Add Ultrapure Water to the microtube of each mix



Ultrapure Water by Contributed by users

#### Step 3.

Add Coral Buffer to each mix.

The Coral Buffer containing 1.5 mM of MgCl2

**■** AMOUNT

 $1 \, \mu l$ : for sample

REAGENTS

1x Coral Buffer by Qiagen

# Step 4.

Add 0.2 mM dNTP

**■** AMOUNT

 $1\,\mu l$ : for sample



# Step 5.

Add control primer (forward and reverse) - HLA-E:

0.04 µl for forward

0.04 µl for reverse

# 5'-3' Sequence:

# **HLA-E** f

CGGGACTGACTAAGGGGCGG

#### HLA-E r

**GTAGCCCTGTGGACCCTCTTAC** 

The size of the fragment is 324 bp

**■** AMOUNT

 $0.04~\mu l: 0.08~\mu M$ 



#### Step 6.

# Add the specific primers:

- MIX 1 (in28-rs11118167 C/ ex29-rs17047660 A; ex26-rs3737002 C/ in26-rs11118131 C)
- MIX 2 (in28-rs11118167 T/ ex29-rs17047660 A; ex26-rs3737002 T/ in26-rs11118131 T)
- MIX 3 (in28-rs11118167 T/ ex29-rs17047660 G; ex26-rs3737002 C/ in26-rs11118131 T)
- MIX 4 (in28-rs11118167\_C/ ex29-rs17047660\_G; ex26-rs3737002\_T/ in26-rs11118131\_C)

0.5µl of each primer (F and R)

We used 0.5  $\mu$ M of SSPs for rs3737002 and rs11118131, 0.6  $\mu$ M for rs11118167 and rs17047660

#### 5'-3' Sequence - rs3737002

CCATTTGCCAGTCCTA**C** 

CCATTTGCCAGTCCTA**T** 

#### 5'-3' Sequence - rs11118131 (reverse primer)

CAAGAAGAAGGGGTGAT**G** 

CAAGAAGAAGGGGTGAT**A** 

#### 5'-3' Sequence - rs11118167

GCCAATATGTGAATATTATTATCTTA**T** 

GCCAATATGTGAATATTATTATCTTA**C** 

# 5'-3' Sequence - rs17047660 (reverse primer)

TTCTGGAGCTGTGCATT**T** 

TTCTGGAGCTGTGCATT**C** 

**■** AMOUNT

0.5 μl: for sample

REAGENTS

✓ Specific Primers Forward and Reverse by Contributed by users

#### Step 7.

Add an aliquot of DNA sample to each tube.

Kepp the plate cold.

**■** AMOUNT

 $1 \mu l$ :

REAGENTS

✓ DNA sample by Contributed by users.

#### Step 8.

Add the Taq DNA polymerase to the mix.

Briefly vortex and centrifuge for some seconds.

AMOUNT

0.056 μl : 0.2 Units/ul



Taq Platinum DNA Polymerase by Invitrogen - Thermo Fisher

# Step 9.

Add an aliquot of the mix to each tube.

Remember, each sample shall be tested in all four mix reactions 1, 2, 3 and 4 (4 tubes)

The final volume for each sample is 10  $\mu$ l.

AMOUNT

9 μl:

# Step 10.

Close well the tubes and centrifuge for some seconds

Put in the thermocycler

# Step 11.

Use:

Thermal cycling begins with  $94^{\circ}\text{C}$  for 5 min and 20 s; followed by 35 cycles, where each cycle began with  $94^{\circ}\text{C}$  for 20 s and ended with  $72^{\circ}\text{C}$  for 40 s. The annealing temperatures were  $59^{\circ}\text{C}$  for the initial 8 cycles,  $57^{\circ}\text{C}$  for the following 7 cycles,  $55^{\circ}\text{C}$  for another 10 cycles and  $53^{\circ}\text{C}$  for the last 10 cycles.

# Warnings

Do not leave the Taq polymerase too long outside the freezer