

Protocol for Bacterial Heterotrophic Production

Jed Fuhrman

Abstract

Citation: Jed Fuhrman Protocol for Bacterial Heterotrophic Production. **protocols.io**

dx.doi.org/10.17504/protocols.io.dmd425


Published: 14 Aug 2015

Guidelines

Materials:

- 10 place filter manifold
- Towers (chilled overnight @ 4C)
- 5% (w/v) Trichloroacetic acid (TCA) (chilled overnight @ 4C) Usually 1L is enough for a single time series set
- 32 15ml centrifuge tubes
- Pump
- Scintillation vials and fluid (UltimaGold)
- Lab supplies

Map:



Protocol

Procedure

Step 1.

Label in duplicate, 4 15ml centrifuge tubes for each depth (3 samples and 1 control)

🔌 NOTES

Jed Fuhrman 10 Aug 2015

One set for Thymidine and another for Leucine.

Procedure

Step 2.

Rinse each tube 3x with 5ml of sample seawater

Procedure

Step 3.

Pour 10ml of sample into rinsed centrifuge tube

Procedure

Step 4.

Kill control samples by adding 0.1ml of filtered formalin

Procedure

Step 5.

Let stand 10 minutes

🕒 DURATION

00:10:00

Procedure

Step 6.

Pipette 5nM of tritiated Thymine or 5nM of tritiated Leucine into sample

🔌 NOTES

Jed Fuhrman 10 Aug 2015

Typically vol. TdR < vol Leu.

Procedure

Step 7.

Dispose of tips into labeled plastic whirlpak bags

Procedure

Step 8.

Incubate samples at seawater temperature for 1 hour

🕒 DURATION

01:00:00

Procedure

Step 9.

During incubation, label scintillation vials

Procedure

Step 10.

10 minutes before incubation time ends set-up Hoeffler box with 0.45µm HA Millipore filters

Procedure

Step 11.

Place chilled towers

Procedure

Step 12.

Following incubation pour samples into respective wells and filter through (<0.3 atm)

📌 NOTES

Jed Fuhrman 10 Aug 2015

Make sure no leaks!!

Procedure

Step 13.

Turn off pump and close wells

Procedure

Step 14.

Add 2ml of 5% TCA to each well

Procedure

Step 15.

Let stand for 2 minutes

🕒 DURATION

00:02:00

Procedure

Step 16.

After 2 minutes, open wells and filter through TCA

Procedure

Step 17.

Rinse towers 3x with 5% TCA

Procedure

Step 18.

Remove towers (pump is still on!)

Procedure

Step 19.

Rinse filter 3x with 5% TCA, enough to cover filters completely

Procedure

Step 20.

Turn off pump

Procedure

Step 21.

Remove filters and place in labelled scintillation vials

Procedure

Step 22.

Add 5ml of UltimaGold scintillation fluid

Procedure

Step 23.

Tightly cap scint vial and shake HARD; make sure the filter ends up in the liquid

Procedure

Step 24.

Let sit for at least 3 hours before final counts, preferably overnight

🕒 DURATION

03:00:00

Procedure

Step 25.

Perform wipe test of surrounding area (see map in guidelines)

Procedure

Step 26.

Read samples in scintillation counter

📌 NOTES

Jed Fuhrman 10 Aug 2015

Currently using program #1 (1 minute read).

Clean up

Step 27.

Rinse towers with diH_2O

Clean up

Step 28.

Rinse top of Hoeffler box with diH_2O

Clean up

Step 29.

Pour radioactive filtrate out of Hoeffler box into waste container carefully

📌 NOTES

Jed Fuhrman 10 Aug 2015

Use a funnel!

Clean up

Step 30.

Throw dry waste into radioactive waste barrel

📌 NOTES

Jed Fuhrman 10 Aug 2015

Not the regular trash!

Clean up

Step 31.

Record isotope usage online