

RNAlater Recipe

Rex Malmstrom

Abstract

This is a storage solution that stabilizes and protects RNA while inactivating RNase

Citation: Rex Malmstrom RNAlater Recipe. [protocols.io](https://doi.org/10.17504/protocols.io.c56y9d)

[dx.doi.org/10.17504/protocols.io.c56y9d](https://doi.org/10.17504/protocols.io.c56y9d)

Published: 12 Oct 2015

Guidelines

Note: If you read the patent application, they refer to an ammonium sulfate concentration of 700g per L. What they really mean is add 700g to 1L (including EDTA and sodium citrate at conc. listed above), so you end up with more than 1L total volume.

Note: The pH should be ~7 using this recipe. The patent states that pH can range from 4-8; the preferred commercial pH is 5.2; I use pH 7 for bacterioplankton samples, but pH can be lowered by adding H₂SO₄.

Note: Since this is a storage solution for RNA, use high quality reagents ("molecular biology" grade) when preparing the buffer.

Protocol

Preparation

Step 1.

Prepare 0.5 M EDTA



. [0.5M EDTA](#)

CONTACT: [Bonnie Poulos](#)

Step 1.1.

Add 186.1 g EDTA disodium salt dihydrate (Mallinckrodt 4931) to 700 ml MilliQ water

Step 1.2.

Stir to dissolve

Step 1.3.

Bring pH to 8.0 with NaOH (50ml of 10M NaOH or 23g NaOH pellets)

Step 1.4.

Bring volume to 1L with MilliQ water

Step 1.5.

Autoclave

Preparation

Step 2.

Prepare 1 M Sodium citrate



PROTOCOL

. [1M Sodium citrate](#)

CONTACT: [Bonnie Poulos](#)

Step 2.1.

Add 294 g Na Citrate dihydrate (Mallinckrodt 0754-12) to 700ml MilliQ water

Step 2.2.

Stir to dissolve

Step 2.3.

Bring volume to 1L with MillQ water

Step 2.4.

Autoclave

RNAlater recipe

Step 3.

In a 2 or 3 L flask mix 1,400 mL MilliQ water with 60 mL 0.5M EDTA and 37.5 mL 1M sodium citrate with stirring

RNAlater recipe

Step 4.

While continuing to stir the solution, add 1.05 kg (= 1,050 g) Ammonium sulfate in 100 g amounts.



REAGENTS

✓ Ammonium sulfate [97061-184](#) by Contributed by users

RNAlater recipe

Step 5.

Stir 30 min on hot plate set to low



DURATION

00:30:00



NOTES

Bonnie Poulos 22 Jun 2015

It is possible that a tiny amount ammonium sulfate will not go into solution. This is okay; the most important thing is to have a saturated ammonium sulfate solution. Crystals should dissolve in ~30min.

RNAlater recipe

Step 6.

Allow solution to cool at room temperature for at least 30min, but overnight on bench is best



DURATION

00:30:00

RNAlater recipe

Step 7.

Filter solution through 0.2 µm filter to remove particles

RNAlater recipe

Step 8.

Store at room temp or at 4°C