

Aug 12,
2019

New Iron Extracting Method from Cattle's Blood for Iron Concentration Analysis V.2

Khairil Asnan¹, Mohd. Sadzan Asram bin Aliming², Aprialdy Idrus², Diana Eka Pratiwi²

¹Biology Department, Universitas Negeri Makassar, ²Chemistry Department, Universitas Negeri Makassar

1 Works for me [dx.doi.org/10.17504/protocols.io.6f2hbqe](https://doi.org/10.17504/protocols.io.6f2hbqe)



ABSTRACT

Cattle's blood component is compressed using centrifuge (5000 rpm for 10 minutes), then supernatant (erythrocytes) collected. Erythrocytes with a ratio of 3:1 (serum:supernatant), treated with initial mixing with NaOH (0.5 M) with an initial ratio 1:1 (v/v) and let it sit for 30 seconds then mixed it by centrifugation about 10 rpm for 30 seconds. Sample then treated by mixing oleic acid (2:1 v/v). The sample then dehydrated by heat about 121°C for a week. Calculated total iron mass was 240000 µg/100 ml whole blood sample (about 14.40% content of the whole sample). Sample purified by furnace using high temperature about 800°C for 2 hours and increased the iron concentration up to 46.30% (m/m%).

MATERIALS

NAME ▾[Sodium Hydroxide](#)**CATALOG #** ▾[View](#)**VENDOR** ▾[Sigma Aldrich](#)

STEPS MATERIALS




NAME ▾[Sodium Hydroxide](#)**CATALOG #** ▾[View](#)**VENDOR** ▾[Sigma Aldrich](#)

MATERIALS TEXT

we also used oleic acid as chelating agent

Sample preparation

1

Compress the whole blood sample from cattle.  6 ml  37 °C  00:10:00



Mini-centrifuge
Centrifuge

Fisher S67601B [↗](#)

Any standard mini centrifuge with adapters for different tube sizes will suffice



 5000 rpm

2

Mix  2 ml blood sample sample with  2 ml NaOH  37 °C  00:00:30

[M] 0.5 Molarity (M)







Sodium Hydroxide
by Sigma Aldrich
[View](#)



dark green solution with strong odor

 10 rpm after we rest the mixed solution for 30 seconds

- 3 Chelate reaction by adding  **4 ml Oleic acid** into the previous mixed sample solution then let it sit for  **00:00:30** . After that, mixed them by  **10 rpm** for  **00:00:30** .



Centrifuge

Benchtop Centrifuge

Eppendorf 5405000441 

Any benchtop centrifuge will suffice



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited