

Running the Thermo Scientific Arktik Thermal Cycler

Dr. Steven Wilhelm, Alyssa Alsante

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Citation: Dr. Steven Wilhelm, Alyssa Alsante Running the Thermo Scientific Arktik Thermal Cycler. **protocols.io**
dx.doi.org/10.17504/protocols.io.in5cdg6

Published: 14 Jul 2017

Protocol

Creating a new protocol

Step 1.

Press NEW from the main screen

NOTES

Alyssa Alsante 27 Jun 2017

This will show a DEFAULT protocol, which is a graphical representation of a typical PCR temperature cycling routine.

Creating a new protocol

Step 2.

Toggle between the different steps with the cursor keys

Creating a new protocol

Step 3.

Press ADD (F1) to add a temperature step to the left of where you want to add it

Creating a new protocol

Step 4.

Press DELETE (F2) to delete a highlighted temperature step

Creating a new protocol

Step 5.

To add advanced options (temperature increment, extend time or slow ramp rate, and temperature gradient) press OPTION (F3) and follow the instructions on the screen.

📌 NOTES

Alyssa Alsante 27 Jun 2017

INC (F1) = add or subtract a temperature increment

EXT (F2) = add or subtract a time extension

RAMP (F3) = slows the ramp rate (input range = 0.1 - 0.4 degrees C/sec.

GRAD (F4) = sets a gradient

Creating a new protocol

Step 6.

Press BACK (F4) to exit without saving

Creating a new protocol

Step 7.

Press SAVE (F5) to save the protocol

Running a protocol

Step 8.

Highlight a folder and press OPEN (F1)

Running a protocol

Step 9.

Highlight the protocol and press OPEN again

Running a protocol

Step 10.

Press START (F1)

Running a protocol

Step 11.

Press OK (F5)

Running a protocol

Step 12.

Specify the SAMPLE (F4) or BLOCK (F5) mode

NOTES

Alyssa Alsante 27 Jun 2017

Sample mode is recommended for most protocols

Running a protocol

Step 13.

Running a protocol

Step 14.

Specify the vessel type by selecting UTW (F4) or COMMON (F5)

Running a protocol

Step 15.

Specify the sample volume

Running a protocol

Step 16.

Press YES (F4) to run the protocol