

# L1 Trace Element Solution

Christa Smith

## Abstract

L1 trace element solution, modified from Bigelow: NCMAM Center. Components to prepare 1 L of trace element solution.

**Citation:** Christa Smith L1 Trace Element Solution. **protocols.io**

[dx.doi.org/10.17504/protocols.io.jvqcn5w](https://dx.doi.org/10.17504/protocols.io.jvqcn5w)

**Published:** 28 Jun 2018

## Guidelines

First, prepare primary stock solutions. Glassware should be acid-washed and rinsed well with dH<sub>2</sub>O before use. If unable to measure small quantities accurately, either scale up the recipe or prepare stock solutions of the individual compounds.

## Before start

Prepare primary stock solutions (in dH<sub>2</sub>O) of the following:

- MnCl<sub>2</sub>·4H<sub>2</sub>O: 178.10 g/L
- ZnSO<sub>4</sub>·7H<sub>2</sub>O: 23.00 g/L
- CoCl<sub>2</sub>·6H<sub>2</sub>O: 11.90 g/L
- CuSO<sub>4</sub>·5H<sub>2</sub>O: 2.50 g/L
- Na<sub>2</sub>MoO<sub>4</sub>·2H<sub>2</sub>O: 19.9 g/L
- H<sub>2</sub>SeO<sub>3</sub>: 1.29 g/L
- NiSO<sub>4</sub>·6H<sub>2</sub>O: 2.63 g/L
- Na<sub>3</sub>VO<sub>4</sub>: 1.84 g/L
- K<sub>2</sub>CrO: 1.94 g/L

Filter sterilize individual solutions through a 0.2 micron PES membrane.

Store frozen in 1 mL aliquots.

## Protocol

### Step 1.

Measure 950 mL dH<sub>2</sub>O into a 1 L glass bottle.

### Step 2.

Dissolve 4.36 g Na<sub>2</sub>EDTA·2H<sub>2</sub>O into the bottle.

### Step 3.

Dissolve 3.15 g FeCl<sub>3</sub>·6H<sub>2</sub>O into the bottle.

### Step 4.

Add 1 mL primary stock of  $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$  (178.10 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 5.**

Add 1 mL primary stock of  $\text{NiSO}_4 \cdot 6\text{H}_2\text{O}$  (2.63 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 6.**

Add 1 mL primary stock of  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$  (23.00 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 7.**

Add 1 mL primary stock of  $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$  (11.90 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 8.**

Add 1 mL primary stock of  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  (2.50 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 9.**

Add 1 mL primary stock of  $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$  (19.90 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 10.**

Add 1 mL primary stock of  $\text{H}_2\text{SeO}_3$  (1.29 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 11.**

Add 1 mL primary stock of  $\text{Na}_3\text{VO}_4$  (1.84 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 12.**

Add 1 mL primary stock of  $\text{K}_2\text{CrO}$  (1.94 g/L in  $\text{dH}_2\text{O}$ ) into the bottle.

---

**Step 13.**

Bring final volume to 1 L with  $\text{dH}_2\text{O}$ .

---

**Step 14.**

Filter sterilize final trace element solution through a 0.2 micron PES membrane.

---

**Step 15.**

Store frozen in 1 mL aliquots.

---