



## Tomato Transformation Media

Version 1

Alex Rajewski<sup>1</sup>

<sup>1</sup>University of California, Riverside

[dx.doi.org/10.17504/protocols.io.x6ffrbn](https://doi.org/10.17504/protocols.io.x6ffrbn)



Alex Rajewski  
University of California, Rive...



### ABSTRACT

This protocol is simply the media recipes for use with my tomato transformation protocol.

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Difco Bacto Agar	156783B	Carolina
Sucrose		
Luria broth powder		
Glycine	GB0235.SIZE.500g	Bio Basic Inc.
Nicotinic acid (NIACIN)	NB0660.SIZE.250g	Bio Basic Inc.
Myo-Inositol	I-525	Gold Biotechnology
Kinetin	K-100	Gold Biotechnology
Pyridoxine HCl (Vitamin B6)	P-780	Gold Biotechnology
Timentin™ Ticarcillin/Clavulanate (15/1)	T-104	Gold Biotechnology
Thiamine HCl	T-260	Gold Biotechnology
trans-Zeatin	Z-105	Gold Biotechnology
Murashige and Skoog with Nitsch Vitamins	MSP29-50LT	Caisson Labs
Murashige & Skoog Basal Salts	MSP01-50LT	Caisson Labs
Agargel	A3301	Millipore Sigma

### MATERIALS TEXT

100mm x 15mm petri dishes (Fisher #FB0875713)

100mm x 20mm petri dishes (VWR #82050-918)

16oz soup containers with lids (Fabri-Kal #9501034 and #9501070)

### BEFORE STARTING

It helps to make up stock solutions of the vitamins, hormones, and antibiotics ahead of time and freeze them.

- Since the concentration of phytohormones and vitamins can vary across medias, we make them both up to a standard concentration of 1mg/mL for easy math.
- Antibiotics are typically used at a single concentration in all of our medias, so they are made up as a 1000x stock. That is, their

actual concentration varies depending on the antibiotic, but in every case you will add 1  $\mu$ L of stock for every 1 mL of media.

- 1 To make 1L, for each media, dissolve the following components to ~800mL of MilliQ water in a 1L beaker

Component:	1/2MSO	MSO-2%	KCMS	2Z	1Z	Rooting
MS Basal Salts	2.15g	4.3g	4.3g	--	--	--
MS Salts with Nitsch Vitamins	--	--	--	4.3g	4.3g	4.3g
Sucrose	10g	20g	30g	20g	20g	30g
myo-Inositol	100mg	100mg	100mg	100mg	100mg	--
KH <sub>2</sub> PO <sub>4</sub>	--	--	200mg	--	--	--
Thiamine HCl (1mg/mL)	2mL	400 $\mu$ L	1.3mL	--	--	--
Pyridoxine HCl (1mg/mL)	500 $\mu$ L	500 $\mu$ L	--	--	--	--
Nicotinic Acid (1mg/mL)	500 $\mu$ L	500 $\mu$ L	--	--	--	--
Glycine	--	2mg	--	--	--	--
2,4-D (1mg/mL)	--	--	200 $\mu$ L	--	--	--
Kinetin (1mg/mL)	--	--	100 $\mu$ L	--	--	--

- 2 Bring volume to 1L and adjust pH ( $\pm$ 0.03) according to the table below:

	1/2MSO	MSO-2%	KCMS	2Z	1Z	Rooting
pH	5.8	5.6	5.5	6.0	6.0	6.0


- 3 Divide the volume into two screw-top jars and add the following gelling agent. Note that not all medias require a gelling agent, and that this mass is the *TOTAL* mass added for 1L. It is not necessary to dissolve this powder; it will dissolve in the autoclave.

	1/2MSO	MSO-2%	KCMS	2Z	1Z	Rooting
Agar	8g	--	--	--	--	--
Agargel	--	--	5.2g	5.2g	5.2g	--
Difco Bacto Agar	--	--	--	--	--	8g

Only half this mass should go into each screw-top jar.

- 4 Autoclave the media, and allow to cool to ~50°. Add any antibiotics and phytohormones (filter-sterile) according to the table below:

	1/2MSO	MSO-2%	KCMS	2Z	1Z	Rooting
Zeatin (1mg/mL)	--	--	--	2mL	1mL	--
Timentin (300mg/mL)	--	--	--	1mL	1mL	--
Antibiotics	No	No	No	Yes	Yes	Maybe

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited