



Feb 27, 2019

Working

Splitting 96 Well Plates for gDNA Extraction and Continuing Culture

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x78frw

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis



Comprehensive Genomic
Editing and Screening
Protocol Updated
02142019.docx

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [Screening Edited iPSC Clones collection](#).

SAFETY WARNINGS









Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

BEFORE STARTING

It takes approximately 1 week for iPSC picked into 96 well plates to be sufficiently confluent for freezing and screening. For screening purposes, a fraction of the cells picked into one well of a 96 well plate will be saved for DNA Extraction and the remaining will be kept in culture or frozen down.

Split cells upon reaching maximum 80% confluence and minimum 40%

- 1 Coat 96 well plate with 50 µl Matrigel per well.
- 2 Incubate at 37 °C for 01:00:00 .
- 3 Prepare plate for expansion by aspirating Matrigel from plate.
- 4 Add 50 µl mTesR1 + 10 uM Rock Inhibitor to appropriate wells.

- 5 Aspirate media from original plate.
- 6 Wash with  200 µl PBS and aspirate.
- 7 Add  25 µl of 0.05% Trypsin.
- 8 Incubate at  37 °C for  00:05:00
- 9 Tap to lift cells from plate.
- 10 Check under microscope to ensure that cells have detached from plate.
- 11 Add  50 µl FBS and tap to mix.
- 12 Transfer  50 µl to 96 well PCR plate, while maintaining the location of each sample (this plate will be used for gDNA extraction).
- 13 Transfer remaining cells (~30uL) to 96 well plate containing mTesR1.
- 14 Incubate at  37 °C .
- 15 After  24:00:00 , complete daily media changes with mTesR1.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited