

Option 2: Bead Cleanup - Alternatively a 2X AMPure XP bead clean up for NEXTflex™ mtDNA-Seq Kit

Nina Orellana

Abstract

Citation: Nina Orellana Option 2: Bead Cleanup - Alternatively a 2X AMPure XP bead clean up for NEXTflex™ mtDNA-Seq Kit. **protocols.io**

dx.doi.org/10.17504/protocols.io.dnq5dv

Published: 14 Sep 2015

Protocol

Step 1.

Add 260 µL of AMPure XP beads to each sample Mix thoroughly.

AMOUNT

250 µl Additional info:

REAGENTS

Agencourt AMPure XP [A63880](#) by [Beckman Coulter](#)

Step 2.

Incubate sample at room temperature for 5 minutes.

DURATION

00:05:00

Step 3.

Place the tube on the magnetic rack at room temperature for 5 minutes or until the supernatant appears clear.

DURATION

00:05:00

Step 4.

Remove and discard clear supernatant taking care not to disturb beads. Some liquid may remain in the tube.

Step 5.

Wash #1: With the tubes on the rack, gently add 500 µL of freshly prepared 80% ethanol to each magnetic bead pellet and incubate at room temperature for 30 seconds. Carefully, remove ethanol by pipette.

DURATION

00:00:30

Step 6.

Wash #2: With the tubes on the rack, gently add 500 µL of freshly prepared 80% ethanol to each magnetic bead pellet and incubate at room temperature for 30 seconds. Carefully, remove ethanol by pipette. **Ensure all ethanol has been removed.**

Step 7.

Remove the tube from the magnetic rack and let dry at room temperature for 3 minutes. Do not overdry the beads.

 DURATION

00:03:00

Step 8.

Resuspend the dried beads with 42 μ L Nuclease-free Water. Mix well by pipetting. Ensure beads are no longer attached to the side of the well.

Step 9.

Incubate resuspended beads at room temperature for 2 minutes.

 DURATION

00:02:00

Step 10.

Place the tube on magnetic rack for 5 minutes or until the sample appears clear.

 DURATION

00:05:00

Step 11.

Gently transfer 40 μ L of clear sample to a fresh microcentrifuge tube.