

# Immunostaining of the isolated lymphoid follicles (ILFs) in the whole small intestine Version 2

Chise Suzuki,Reiji Aoki,Ayako Aoki-Yoshida

## Abstract

The immunostaining detection of isolated lymphoid follicles (ILFs) was performed according to the method of McDonald and Newberry [BioTechniques. 2007; 43(1): 50-56], with some modifications. Instead of mounting the tissue samples on a mounting plate, each segment was treated in a 15 ml tube in each step. DAB stained intestines were finally mounted on a glass slide using an aqueous mounting medium.

**Citation:** Chise Suzuki,Reiji Aoki,Ayako Aoki-Yoshida Immunostaining of the isolated lymphoid follicles (ILFs) in the whole small intestine. **protocols.io**

dx.doi.org/10.17504/protocols.io.jr3cm8n

**Published:** 12 Sep 2017

## Materials

Hanks' Balanced Salt solution, 10 ×, Modified, without calcium, magnesium or sodium bicarbonate [H4385-100ML](#) by [Sigma Aldrich](#)

Avidin/Biotin Blocking Kit [SP-2001](#) by [Vector Laboratories](#)

Purified Rat Anti-Mouse CXCR5 [551961](#) by [BD Biosciences](#)

Anti-Rat Ig HRP detection kit [551013](#) by [BD Biosciences](#)

✓ Mount-Quick Aqueous [View](#) by Contributed by users

## Protocol

### Preparation

#### Step 1.

Open small intestine vertically

#### Step 2.

Rince the intestine with cold PBS in petri dish

#### Step 3.

Divide a small intestine into four segments (A to D, proximal to distal).

#### **Step 4.**

Place each segment into a 15 ml tube filled with 10 ml warmed (37°C) Hank's balanced salt solution (HBSS) containing 5 mM EDTA at 37°C

#### **Removal of the epithelial cells**

#### **Step 5.**

Incubate tubes on shaker at 37°C for 15 min (shaking 90 rpm) and then shake the tubes vigorously 20 times by hand.

#### **Step 6.**

Discard supernatant containing epithelium and filter remaining the intestinal segment through a gauze.

#### **Step 7.**

Put retained the intestinal segment in the 15 ml tube containing 10 ml HBSS containing 5 mM EDTA.

#### **Step 8.**

Repeat steps 5-7 for 3 times.

#### **Fixation**

#### **Step 9.**

Rince the segment with cold PBS and submerge the segment in 10% formalin-buffered saline and place at 4°C at least 1h.

#### **Step 10.**

Wash the segment for 5 min with 50 mM Tris-HCl (pH 7.2) buffered saline (150 mM NaCl) (TBS) containing 0.5% Triton X-100 for 3 times.

#### **Quenching**

#### **Step 11.**

Add 1% H<sub>2</sub>O<sub>2</sub> in methanol and shake 15 min at room temperature (rt).

#### **Step 12.**

Wash for 5 min with TBS containing 0.05% Tween 20 (TBST) for 3 times.

## Blocking

### Step 13.

Add 1 drop Avidin D solution (Avidin/Biotin Blocking Kit) in 250 µl PBS containing 1% BSA and incubate 30 min at rt.

Then wash the segment with TBS.

## Primary antibody (anti-CXCR5)

### Step 14.

Add 1 drop of biotin solution (Avidin/Biotin Blocking Kit) and 2.5 µl of primary antibody\* in 250 µl TBST (\* Rat anti-mouse CXCR5 (clone 2G8, BD Biosciences, San Diego, CA) and incubate overnight at 4°C.

Wash the segments for 5 min with PBS 3 times.

## Secondary antibody (Anti Rat Ig HRP)

### Step 15.

Add 2 µl biotinylated anti-Ig secondary antibody (Anti-Rat Ig HRP detection kit, Pharmingen) in 400 µl antibody diluent and incubate for 30 min at rt.

Wash for 2 min with PBS 3 times.

### Step 16.

Add 350 µl Streptavidin-HRP and incubate for 30min at rt.

Wash for 2 min with PBS 3 times.

## Detection

### Step 17.

Drain PBS. Add 300 µl 3,3'-diaminobenzidine tetrahydrochloride (DAB) solution (1 drop DAB/1ml DAB buffer).

Incubate for 5 min or until the desired color intensity is reached.

Wash for 2 min with PBS 3 times.

Rince the segment with water for 2 min 3 times

## Mounting

### Step 18.

Mount the stained segments of intestine on a slide glass and fixed with Mount-Quick Aqueous.