

3' Azide labelling with TdT

Benjamin Emert

Abstract




Protocol for add azide modified nucleotide to 3' end of unmodified DNA oligonucleotide using terminal transferase. Protocol adapted from [Winz et al. 2015](#).

Citation: Benjamin Emert 3' Azide labelling with TdT. **protocols.io**

dx.doi.org/10.17504/protocols.io.qnmdvc6

Published: 01 Jun 2018

Materials







-  Terminal Transferase - 2,500 units [M0315L](#) by [New England Biolabs](#)
-  ethanol by Contributed by users
- N6-(6-Azido)hexyn-3dATP [NU-1707](#) by [Jena Bioscience](#)
-  3M Sodium Acetate solution by Contributed by users

Protocol

TdT reaction

Step 1.

Set up TdT reaction

-  [AMOUNT](#)
2.5 µl Additional info: 10x Terminal transferase buffer
-  [AMOUNT](#)
10 µl Additional info: 2.5mM CoCl₂
-  [AMOUNT](#)
2 µl Additional info: 400µM template oligo
-  [AMOUNT](#)
1 µl Additional info: 10mM N6(6-azido)hexyl-3'-dATP
-  [AMOUNT](#)
2.5 µl Additional info: recombinant TdT (20U/µL)
-  [AMOUNT](#)
7 µl Additional info: nuclease free water

[NOTES](#)

Benjamin Emert 01 Jun 2018

Probably useful to scale up reaction.

TdT reaction

Step 2.

Incubate at 37°C overnight.

TEMPERATURE

37 °C Additional info: Incubate on thermal cycler

TEMPERATURE

50 °C Additional info: thermal cycler lid temperature

Ethanol precipitate to cleanup

Step 3.

Add 1/10th volume sodium acetate and 3x volume ethanol to precipitate azide labelled oligo.

AMOUNT

2.5 µl Additional info: 3M sodium acetate

AMOUNT

82.5 µl Additional info: 200 proof ethanol

NOTES

Benjamin Emert 01 Jun 2018

If you want to get rid of excess N6-(6-azido)hexyl-3'-dATP, you may want to use 1/2 volume 7.5M ammonium acetate.

Ethanol precipitate to cleanup

Step 4.

Place solution at -80°C for 1 hour to precipitate.

TEMPERATURE

-80 °C Additional info:

NOTES

Benjamin Emert 01 Jun 2018

Fine to incubate at -20°C or -80°C for longer.

Ethanol precipitate to cleanup

Step 5.

Centrifuge for 10 minutes to precipitate

TEMPERATURE

4 °C Additional info:

Ethanol precipitate to cleanup

Step 6.

Remove supernatant. Careful not to disturb pellet.

Ethanol precipitate to cleanup

Step 7.

Wash sample with 500 µL 70% ethanol

AMOUNT

500 µl Additional info: 70% ethanol

Ethanol precipitate to cleanup

Step 8.

Centrifuge for 2-5 minutes.

Ethanol precipitate to cleanup

Step 9.

Remove supernatant. Careful not to disturb pellet.

Ethanol precipitate to cleanup

Step 10.

Centrifuge for 1 minute to collect supernatant from the lid and walls.

Ethanol precipitate to cleanup

Step 11.

Remove supernatant then dry pellet under vacuum or ambient air.

Ethanol precipitate to cleanup

Step 12.

Add 1.8 µL nuclease free water to dissolve dry pellet. Let sit at room temp for ≥20 minutes.

AMOUNT

1.8 µl Additional info: nuclease free water

TEMPERATURE

25 °C Additional info: incubate at room temperature

NOTES

Benjamin Emert 01 Jun 2018

Adding 1.8 µL and assuming ~10% loss. Can be more precise and measure final concentration but beware of residual dATP.

Ethanol precipitate to cleanup

Step 13.

Store azide modified oligo at -20°C indefinitely.

TEMPERATURE

-20 °C Additional info: