



# untitled protocol

#### Yuan Yao

#### **Abstract**

Citation: Yuan Yao untitled protocol. protocols.io

dx.doi.org/10.17504/protocols.io.qaadsae

Published: 22 May 2018

## **Materials**



IPTG 12481 by Gold Biotechnology

#### **Protocol**

#### Step 1.

Incubate BW25113/pRstA-GFP cells in 50 ml LB medium with chloromycetin (100 mg/ml), and expression was induced by adding IPTG (final concentration of 1.0 mM).

#### Step 2.

Harvest by centrifugation, the cells were resuspended in 12 ml Lysis/Equilibration buffer including 120 µl lysozyme (10 mg/ml),

#### Step 3.

Place on ice for 30 min, and then extracted by sonication method.

#### Step 4.

Centrifuge at 4 °C 12000 rpm for 5 min.

## Step 5.

Transfer supernatant to new 1.5 ml epp tube.

## Step 6.

Incubate with His-Select Ni-NTA Agarose (Thermo Fisher Scientific, USA) over night at 4°C.

#### Step 7.

The next day, the supernatant was loaded onto Ni-NTA affinity resin (Qigen, Germany).

## Step 8.

Proteins were examined for purity by SDS-PAGE at 200 V for 45 min. Stain gel with extant blue for 15 min.

# Step 9.

The pure protein were pooled and dialyzed.

# Step 10.

The extraction of proteins was analyzed by BCA protein concentration quantitative method.