

Digestion NotI

Antoine ZILLER

Abstract

Citation: Antoine ZILLER Digestion NotI. [protocols.io](https://doi.org/10.17504/protocols.io/ecubaww)

[dx.doi.org/10.17504/protocols.io/ecubaww](https://doi.org/10.17504/protocols.io/ecubaww)

Published: 05 Jan 2016

Protocol

Step 1.

Add:

- nuclease-free water qsp 50 μ L > 47 μ L
- 5 μ L 10X Buffer O
- 1 μ L DNA (1 μ g/ μ L)
- 1 μ L NotI 1 (10U/ μ L) *

- Mix gently and spin down for a few seconds.

Step 2.

Incubate at 37°C for 1-16 hours. The digestion reaction may be scaled either up or down

> It's possible to process extended incubation by adding 0.25U/ μ g of DNA in 50 μ L of reaction volume.

 **DURATION**

01:30:00

Step 3.

incubation at 80°C for 20 min.

Step 4.

1. Measure the volume of the DNA sample.

> 50 μ L

2. Add 1/10 volume of sodium acetate, pH 5.2, (final concentration of 0.3 M) - These amounts assume that the DNA is in TE only; if DNA is in a solution containing salt, adjust salt accordingly to achieve the correct final concentration.

> 5 μ L

3. Mix well.

4. Add 2 to 2.5 volumes of cold 100% ethanol (calculated after salt addition).

> 110 (2V)

5. Mix well.

6. Place on ice or at -20 degrees C for >20 minutes.

7. Spin a maximum speed in a microfuge 10-15 min.
8. Carefully decant supernatant.
9. Add 1 ml 70% ethanol. Mix. Spin briefly. Carefully decant supernatant.
10. Air dry or briefly vacuum dry pellet.
11. Resuspend pellet in the appropriate volume of TE or water.