



## *E. coli* Optical Quantification

Version 2

Kenneth Schackart<sup>1</sup>, Kattika Kaarj<sup>1</sup>

<sup>1</sup>University of Arizona

[dx.doi.org/10.17504/protocols.io.xu3fnyn](https://doi.org/10.17504/protocols.io.xu3fnyn)

481b Laboratory

 Kenneth Schackart 

### ABSTRACT

This procedure details how to culture *E. coli* K-12 in suspension and quantify cultured bacteria using a spectrometer.

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working




### GUIDELINES

Labcoat and gloves must be worn at all times.

### MATERIALS TEXT

- Gloves
- Lyophilized *E. coli* K-12
- BHI Broth
- 15 mL centrifuge tubes
- Pipette and tips
- 70% Ethanol solution
- Kimwipes®
- Mini Spectrometer
- Cuvettes
- Computer with OceanView software

### Inoculate Culture

- 1 Pipette  BHI broth into 15 mL centrifuge tube.
- 2 Place tube in incubator or water bath to warm.
- 3 Add  lyophilized *E. coli* K-12 to warm BHI broth.
- 4 Close lid and shake gently.
- 5 Incubate at  for several hours.

## Quantify Concentration

- 6 Dispense  1 ml BHI broth into cuvette.

**NOTE**

This should be broth *without* bacteria and will serve as your reference for calculating absorbance.

- 7 Carefully insert cuvette into spectrometer.

**NOTE**

Wipe cuvette with delicate task wipes if necessary.

- 8 Measure intensity at 600 nm.

- 9 Repeat Steps 6-8 with the samples cultured for different durations.

**NOTE**

Avoid pipetting any precipitation at the bottom of the culture tube. Draw only supernatant.

- 10 Using the BHI broth without bacteria as  $I_0$ , calculate the absorbance of each of the samples.

- 11 Plot Absorbance as a function of time cultured. **Include this plot, and the table of intensities and calculated absorbances in the 'Results' section of your lab report.**



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited