

# Zobell Broth

Sarah Hessen-Schmidt

## Abstract

**Citation:** Sarah Hessen-Schmidt Zobell Broth. [protocols.io](https://protocols.io)

[dx.doi.org/10.17504/protocols.io.e7wbhpe](https://dx.doi.org/10.17504/protocols.io.e7wbhpe)

**Published:** 21 Jun 2016

## Protocol

### Step 1.

Add 40.25g of Zobell Marine Broth to 1000ml MilliQ. Mix thoroughly. Heat if necessary to dissolve.

### Step 2.

Place autoclave tape on container and autoclave. After autoclaving, store in fridge.