

# Developing a Titration Kill Curve (G418, Hygromycin B and Puromycin)

Sean Seaver

## Abstract

Developing a Titration Kill Curve (G418, Hygromycin B and Puromycin)

**Citation:** Sean Seaver Developing a Titration Kill Curve (G418, Hygromycin B and Puromycin). **protocols.io**

dx.doi.org/10.17504/protocols.io.ci9uh5

**Published:** 03 Nov 2014

## Before start

Each mammalian cell line has a different sensitivity. Before experimentation, you should determine the optimal concentration of your antibiotic by developing the Titration Kill Curve for your cells.

## Materials

- G418 Disulfate [LGB-418](#) by [P212121](#)
- Hygromycin B [GB-T005](#) by [P212121](#)
- Puromycin Dihydrochloride [RP-P33020](#) by [P212121](#)

## Protocol

### Step 1.

Split non-transduced, confluent cells 1:5 in 10 ml DMEM + 10% FBS media.

 DURATION

00:05:00

### Step 2.

Transfer 0.5 mL cell suspension into 24-well plate containing 500 µl of (media + drug).

- [G418 Sulfate](#)
- [Hygromycin B](#)
- [Puromycin Dihydrochloride](#)

 DURATION

00:08:00

### Step 3.

Examine viability every 2 days.

 DURATION

48:00:00

#### Step 4.

Culture for 14 days. Replace the media containing antibiotic every 3 days.

#### DURATION

12:00:00

#### Step 5.

Use the lowest concentration of your antibiotic that begins to give massive cell death in 3 days and kills all the cells within two weeks.

A general starting point is usually 400 mg/ml G418 for HeLa cells and 200 mg/ml hygromycin for CHO cells.

In mammalian cells the optimal level of puromycin is typically around 1 mg/ml.

HeLa cells are often selectable with 500 mg/ml G418, 500 mg/ml hygromycin, or 2.5 mg/ml puromycin, and SHSY-5Y cells are often selectable with 600 mg/ml G418 or 200mg/ml hygromycin.

#### REAGENTS

 Hygromycin B [GB-T005](#) by [P212121](#)