

Culturing Euplotes crassus to high densities using E. coli as the only food source Version 10

Rachele Cesaroni

Abstract

Citation: Rachele Cesaroni Culturing Euplotes crassus to high densities using E. coli as the only food source. protocols.io

dx.doi.org/10.17504/protocols.io.hicb4aw

Published: 31 Mar 2017

Protocol

Step 1.

Grow 1 L culture of E. coli in Luria broth to saturation overnight (usually we dilute previous bacterial culture 1:100).

📌 NOTES

Rachele Cesaroni 31 Mar 2017

We have been using strain HT115, but any strain of E. coli will likely do. Do not use antibiotics.

Step 2.

For 1 L Euplotes crassus culture, pellet E. coli from 200 ml of culture (4000 rpm for 10 minutes). The remaining bacteria can be stored in their 1 L flask at 4 °C for at least a month and used to feed Euplotes crassus as necessary.

Step 3.

Wash the pellet once with ddH₂O and pellet it again at 4000 rpm for 10 minutes.

Step 4.

After discarding most of the excess water, resuspend the bacteria (e.g. with a micropipette using a 1 ml tip) before adding them to the Euplotes crassus culture.

Step 5.

Euplotes crassus cells typically consume all the bacteria after 2 or 3 days at 24 °C with aeration system, reaching a density of 3000 cells/ml.