

Cas9 RNP nucleofection for cell lines using Lonza 4D Nucleofector

Mark Dewitt & Julia Wong

Abstract

Citation: Mark Dewitt & Julia Wong Cas9 RNP nucleofection for cell lines using Lonza 4D Nucleofector. **protocols.io**
dx.doi.org/10.17504/protocols.io.dm649d

Published: 13 Aug 2015

Protocol

Prepare RNP mix

Step 1.

Bring 100 pmol of Cas9 to a final volume of 5 μ L using Cas9 buffer (20 mM HEPES-KOH pH 7.5, 150 mM KCl, 10% glycerol, 1 mM TCEP). For 40 μ M stock: 2.5 μ L.

■ ANNOTATIONS

Jacob Corn 02 Mar 2017

Preparing RNP mix notes: Cas9-NLS is stored in -80°C, sgRNAs are prepped by runoff transcription, Cas9 buffer is kept in the TC hood and must be kept sterile.

Prepare RNP mix

Step 2.

Bring 120 pmol sgRNA to a final volume of 5 μ L using Cas9 buffer. This means you will need a minimum sgRNA concentration of 24 μ M.

Prepare RNP mix

Step 3.

Add Cas9 to sgRNA slowly while swirling pipette tip, should take 30s to 1 minute.

Prepare RNP mix

Step 4.

Allow RNP to form for 10-20 minutes.

🕒 DURATION

00:20:00

Prepare Cells

Step 5.

Count cells. (Trypsinize as needed.)

Prepare Cells

Step 6.

For each nucleofection, pipette 200k cells into a 15 mL conical.

Prepare Cells

Step 7.

Spin 100 x g for 10 minutes to pellet cells softly. While the cells are spinning, prepare plate and

cuvette.

 **DURATION**

00:10:00

Prepare Cells

Step 8.

Prepare a 12-well-plate with 1mL media per well, and pre-warm in the incubator.

Nucleofection

Step 9.

Prepare and label wells on 20uL nucleofection strips. Configure Lonza 4d using recommended cell-type program.

Nucleofection

Step 10.

Pipette off media from cells, gently but completely, using a P200. The pellet is very soft so be careful.

Nucleofection

Step 11.

Resuspend cells in 20 μ L of nucleofector solution (usually SF media) using a P200.

Nucleofection

Step 12.

Add the entire 10 μ L RNP mix to the 20 μ L resuspension and mix.

Nucleofection

Step 13.

Add 1uL of 100uM donor DNA (100 pmoles) and mix well.

Nucleofection

Step 14.

Add nucleofection mixes to the multiwell cuvette, and cap. Pay attention to the orientation of the cap and cuvette in the nucleofector, which is noted in the manufacturer's instructions.

Nucleofection

Step 15.

Insert cuvette into nucleofector and zap.

Nucleofection

Step 16.

Allow cells to sit in nucleofection strips for 10 minutes post-nucleofection. This is supposed to increase efficiency.

 **DURATION**

00:10:00

Nucleofection

Step 17.

Add 80uL of pre-warmed media to each well.

Nucleofection


Step 18.

Pipette mixture out with a P200 into your pre-warmed 12-well plate. This should get the vast majority of cells, but if you wish, you may wash out the rest with media from the same well, chemistry-style.

Nucleofection

Step 19.

Allow cells 24 hours to settle and recover before attempted downstream analysis. Consider including un-zapped controls to test viability.

 DURATION
24:00:00