

# Efficiency of Bacteriophage Plaquing (EOP) Assay Version 2

Jessica Sacher

## Abstract

Efficiency of plaquing (EOP) assay for *C. jejuni* phages

Jessica Sacher (with help from Stephanie Guerra)

**Citation:** Jessica Sacher Efficiency of Bacteriophage Plaquing (EOP) Assay. **protocols.io**

dx.doi.org/10.17504/protocols.io.mahc2b6

**Published:** 18 Dec 2017

## Protocol

### Step 1.

For each strain, streak jejuni cells from frozen onto NZCYM 1.5% agar plates, grow 2 days in microaerobic incubator.

### Step 2.

Restreak cells onto a fresh NZCYM plate and grow 1 day.

### Step 3.

Harvest cells by flooding plate with 5 mL NZCYM broth, mixing with a spreader and collecting resultant suspension.

### AMOUNT

5 ml Additional info: NZCYM broth

### Step 4.

Mix well by inversion and set OD<sub>600</sub> for each strain to 0.35 (or to 0.0035 – this has also worked well) in 5 mL NZCYM broth in an empty Petri dish (regular size).

### Step 5.

Incubate at 37°C for 4 h at 100 rpm (place Petri dishes in a large glass beaker inside of an autoclave bag and fill with trigas mixture – beaker should also contain 50 mL water in a flask which is taped to the side of the beaker, which helps to keep atmosphere moist).

### Step 6.

During 4-h incubation:

- Melt NZCYM 0.6% agar (microwave) and keep at 55°C until ready to harvest cells (let sit at least 30 mins at 55C prior to use in order to ensure it is not too hot).
- Serially dilute phages (10-fold dilutions, usually down to 10<sup>-5</sup> or 10<sup>-6</sup>) and store on ice until

needed.

- Label (including spaces to spot dilutions: 0, -1, -2, -3, -4, -5) and pre-dry (open inside biosafety cabinet for 30 min) NZCYM plates (one per phage/strain combination).

### **Step 7.**

After the 4-h incubation, mix cultures by swirling and measure OD<sub>600</sub>. Set OD of each strain to 0.3 and add 250 µL of this suspension to a 15-mL Falcon tube (0.5 OD x 200 µL has also worked well).

### **Step 8.**

Add 5 mL molten agar to each tube, mix thoroughly but gently (avoid bubbles) by inversion. Immediately pour the entire contents of each tube onto an NZCYM plate.

### **Step 9.**

Let solidify 10-20 mins.

### **Step 10.**

Spot dilutions of phages (10 µL x 2 spots per dilution) as well as undiluted phage stock and let plates dry with lids off until spots are absorbed (15-20 min).

### **Step 11.**

Incubate plates upside down in microaerobic incubator for 18-24 h. Check plates for plaques and image using an iPhone (or appropriate substitute). Count isolated plaques and record.

### **Step 12.**

To convert plaque numbers to PFU/mL, multiply the number of isolated plaques by the inverse of the dilution factor (ie. If 5 plaques are recorded on the -4 dilution spot, multiply 5 by 10<sup>4</sup>, then multiply this number by 100 to account for pipetting only 10 µL – this would give 5 x 10<sup>6</sup> PFU/mL). Shortcut: number of plaques (5, in this example) x 10<sup>(dilution number plus 2)</sup> = 5 x 10<sup>(4+2)</sup> = 5.0 x 10<sup>6</sup> PFU/mL.

### **Step 13.**

Graph PFU/mL vs. jejuni strain to compare plaquing efficiency of a phage on multiple strains. Alternatively, convert each PFU/mL value to a percentage of a reference (e.g. wild type) strain.