

# Protein renaturation from pellet samples

Snehadri Sinha

## Abstract

After cell lysis and centrifugation the possibly misfolded proteins can be refolded with this protocol.

**Citation:** Snehadri Sinha Protein renaturation from pellet samples. **protocols.io**

dx.doi.org/10.17504/protocols.io.gzkbx4w

**Published:** 17 Jan 2017

## Protocol

### Step 1.

Resuspend the samples in 50 mM Tris-HCl at pH 7.4, 4 M urea

### Step 2.

Sonicate samples for four minutes, 5 sec on, 5 sec off, at 30 % Amp to wash

### Step 3.

Centrifuge 12 000 g for 10 min, RT

### Step 4.

Collect supernatant

### Step 5.

Repeat the washing three times, collect supernatant

### Step 6.

Incubate the supernatant at 37 °C for minimum 1 hour

### Step 7.

If your sample is cloudy after this add 8 M urea until it becomes clear, incubate for 1 h at 37 °C

### Step 8.

Dialyze your denaturated protein samples first at a 1:10 ratio against buffer containing 10 % v/v glycerol and 0.1 mM EDTA (refolding buffer), 20 mM HEPES pH 7.4 at 4 °C for 4 h

### Step 9.

Optional: Concentrate your sample

### Step 10.

Second dialysis with a 1:100 ratio against the same buffer at 4 °C for 16 h or O/N

### Step 11.

Concentrate samples for further analysis