



Oct 18, 2019

Double Digestion and Dephosphorylation of Plasmid

iGEM Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.79phr5n](https://doi.org/10.17504/protocols.io.79phr5n)

iGEM Dusseldorf ⚡

1

Mix the following components gently:

Components	ng/μL
Plasmid DNA Tube X	
Component	Volume (μL)
Sterile MilliQ Water	Fill up to 30.0 μl
10x FastDigest Buffer	3.0
Plasmid DNA	Total ng of DNA/Plasmid DNA concentration ng/μl
FastDigest Enzyme 1	1.0
FastDigest Enzyme 2	1.0
FastAP Alkaline Phosphatase	1.0
Total Volume	30.0
Total ng of DNA	1000.0

2 Incubate at enzyme suitable conditions ⌚ 00:30:00 - ⌚ 01:00:00 and temperature



Find out suitable conditions on website of particular science company

3 Inactivate the enzymes at suitable conditions and temperature



Find out suitable conditions on website of particular science company

4 Run gel at 80-150 V until the dye line is approximately 75-80% of the way down the gel. A typical run time is about 1-1.5 hours, depending on the gel concentration and voltage.

5 Cut out digested band (~150 mg) and transfer to pre-weighed tube. Check weight.

6 Purify the band



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited