

Dephosphorylation of 5'-ends of DNA using CIP (M0290)

Version 2

New England Biolabs

Abstract

Protocol for Dephosphorylation of 5'-ends of DNA using CIP in Restriction Enzyme Reaction. Uses the Calf Intestinal Alkaline Phosphatase (CIP - M0290)

Citation: New England Biolabs Dephosphorylation of 5'-ends of DNA using CIP (M0290). **protocols.io**

dx.doi.org/10.17504/protocols.io.iqwcdxe

Published: 08 Aug 2017

Guidelines

Dephosphorylation of 5'-ends of DNA in Restriction Enzyme Reaction

- The phosphate can be added directly into the digestion reaction during or after DNA digestion
- CIP is active in all NEB restriction enzyme buffers
- DNA purification is required before ligation

Materials

🐼 Alkaline Phosphatase, Calf Intestinal (CIP) [M0290](#) by [New England Biolabs](#)

Protocol

Step 1.

Prepare a 20 µl reaction as follows:

DNA	1 pmol of DNA ends*
CutSmart® Buffer (10X)	2 µl
CIP	1 unit
H ₂ O, purified	to 20 µl**

📌 NOTES

New England Biolabs 30 Jun 2017

* *Note: 1 pmol of DNA ends is about 1 µg of a 3 kb plasmid.*

*** Scale larger reaction volumes proportionally.*

Step 2.

Incubate at 37°C for 30 minutes.

 **DURATION**

01:00:00

Step 3.

Purify DNA by gel purification, spin-column or phenol extraction.