



MALDI-TOF MS preparation for identification of mosquitoes

Version 3

Andrea L. Lawrence¹, Maureen Laroche², Jana Batovska³, Cameron E. Webb¹, Stacey E. Lynch³, Mark J. Blacket³, Jan Šlapeta⁴, Philippe Parola²

¹Marie Bashir Institute of Infectious Diseases and Biosecurity, University of Sydney, Sydney, New South Wales 2006, Australia; Medical Entomology, NSW Health Pathology, ICPMR, Westmead Hospital, Westmead, New South Wales 2145, Australia, ²Aix Marseille University, IRD, AP-HM, SSA, VITROME, IHU-Méditerranée Infection, 19-21 Boulevard Jean Moulin 13005 Marseille, France, ³Agriculture Victoria, AgriBio, Centre for AgriBioscience, Bundoora, Victoria 3083, Australia, ⁴Sydney School of Veterinary Science, Faculty of Science, University of Sydney, Sydney, New South Wales 2006, Australia

dx.doi.org/10.17504/protocols.io.vnee5be



ABSTRACT

MALDI-TOF MS specimen preparation for rapid identification of Australian mosquitoes

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

There are 3 main steps involved with preparing and launching a mosquito sample for MALDI-TOF analysis:

- Preparing the matrix

The matrix solution ionises the proteins in the sample and is applied to the plate after the sample. The matrix must be very acidic.




- Dissecting and crushing the sample
- Preparing and launching the plate

MATERIALS








NAME ▾	CATALOG # ▾	VENDOR ▾
1.5 mL Eppendorf tubes		
Acetonitrile	AC1400.SIZE.1L	Bio Basic Inc.
Trifluoroacetic acid (TFA)	TC8960.SIZE.100mL	Bio Basic Inc.
Water, uHPLC grade		
α-cyano-4-hydroxycinnamic acid	C8982 SIGMA	Sigma Aldrich
Formic acid, 70%		
Glass beads, acid washed ≤106µm	G4649-500G	Sigma Aldrich

Preparing the matrix

- Using a small spatula, put approx. **50 µl matrix powder (α-cyano-4-hydroxycinnamic acid)** in the Eppendorf tube and return the stock powder to the fridge immediately.
- Add **500 µl 100% ACN (acetonitrile)**, **475 µl HPLC grade water** and **25 µl 100% trifluoroacetic acid** to the Eppendorf tube.

- 3 Put Eppendorf tube with solution in an ultrasonic bath for  **00:10:00** to homogenise the solution (yellow matrix powder will still be visible); then centrifuge the solution for  **00:10:00** at 13 000 rpm/ 20 784 x g.
- 4 Retrieve the supernatant using a pipette, taking care to avoid taking up any powder. Discard the powder. The matrix can be kept at  **4 °C** in the dark for up to 2-3 days when creating high quality spectra for reference database. It can be kept for up to a week when using it for simple identification when high quality spectra are not vital.

Preparing the sample

- 5 If using live mosquitoes, put the mosquitoes in the freezer ( **-20 °C**) for at least  **00:01:00** prior to dissection.
- 6 Using two sets of forceps, carefully remove the legs from the mosquito, using one pair to steady the specimen and the other to remove the legs. Place the legs in a labelled and sterile Eppendorf tube. The legs will stick to the sides of the tube due to static so try to put the legs at the bottom of the tube.
- 7 If more than one species are being processed, decontaminate the forceps with 70% ethanol between species.
- 8 Spin down the legs for  **00:01:00** at 13 000 rpm/ 20 784 x g.
- 9 Add  **15 µl 50% CAN** and  **15 µl 70% formic acid** to each tube then add glass beads (<106µm, acid washed). Estimate the volume, making the amount of glass beads just less than half of the volume of liquid.
- 10 Place tubes in the TissueLyser II (Qiagen) canisters, making sure each side is balanced. NB: Both the tubes and the canisters need to be balanced. Make sure the top and bottom lids correlate and that the heavy sides are both either facing inward or both outward.
- 11 Crush the legs using a cycle of 3 x  **00:01:00** at 30 Hz, letting the sample rest for approx.  **00:00:15** between each minute.
NB: the legs will not be fully crushed in the tube; too much protein will inhibit the reaction as in PCR.

Preparing the plate



- 12 Load  **1 µl of each sample** on the plate in quadruplicate keeping within the spots.
- 13 Allow the sample to evaporate then load  **1 µl matrix solution** over the sample, allow to dry completely.
- 14 Place the plate in a cool, dark (photosensitive matrix) place until it's time to launch it. NB: the prepared plates can be left for max. 2-3 days in a dark place at room temperature before affecting the spectra. If using for simple identification – it can be kept for a longer period.

Plate is ready to be launched in MALDI-TOF MS machine.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited