

# Adsorbing Viruses on TEM Grids

Jennifer Brum

## Abstract

**Purpose:** This protocol describes how to adsorb viruses onto TEM (transmission electron microscopy) grids. The sample is allowed to sit on a hydrophilic grid and viruses adsorb onto the surface of the grid. This technique is generally used for viral lysates with high concentrations of viruses. For natural samples, use the protocol "[Quantitatively Depositing Viruses onto TEM Grids using an Airfuge](#)".

Note: If purification of the viruses is necessary, refer to Ackermann and Heldal (2010) for various options.

**Citation:** Jennifer Brum Adsorbing Viruses on TEM Grids. **protocols.io**  
dx.doi.org/10.17504/protocols.io.dar2d5

**Published:** 21 Jan 2016

## Guidelines

### Materials required:

- glow discharge apparatus
- EM grid-grade tweezers (2 or more pairs)
- EM grids (formvar coated 200 mesh copper)

**Note:** If purification of the viruses is necessary, refer to Ackermann and Heldal (2010) for various options.


## Protocol

### Step 1.

On the same day as grid prep, use glow discharge to render grids hydrophilic.

### Step 2.

Holding the grid in the tweezers, place 5 µl of viral lysate onto the shiny side of the grid (Figure 1).



### Step 3.

Let the sample sit for 3 minutes to allow the viruses to adsorb to the grid.

 DURATION

00:03:00

**Step 4.**

Stain the viruses using the protocol “[Positive and Negative Staining of Viruses on TEM Grids](#)”.