

Oct 23, 2019

Lentiviral transduction of iPSCs with sgRNAs and sgRNA libraries

Ruilin Tian¹, Jason Hong¹, [Sydney Sattler¹](#), [Martin Kampmann¹](#)¹University of California, San Francisco[1](#) *Works for me* dx.doi.org/10.17504/protocols.io.8dfhs3n[Neurodegeneration Method Development Community](#) [KampmannLab](#)**Martin Kampmann**
University of California, San Francisco

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Tian et al (2019). CRISPR Interference-Based Platform for Multimodal Genetic Screens in Human iPSC-Derived Neurons. Neuron pii: S0896-6273(19)30640-3. [Epub ahead of print] PubMed PMID: 31422865.

Lentiviral transduction of
iPSCs with sgRNAs and
sgRNA libraries.pdf

MATERIALS

NAME	CATALOG #	VENDOR
DMEM, high glucose	11965092	Thermo Fisher Scientific
Opti-MEM™ I Reduced Serum Medium	31985070	Thermo Fisher Scientific
TransIT®-Lenti Transfection Reagent	MIR 6600	Mirus Bio
Lentivirus Precipitation Solution	VC125	ALSTEM Cell Advancements

MATERIALS TEXT



12 ml Luer Lock Syringe
Syringe

NORM-JECT® 4100.X00V0 [↗](#)



Filter, 0,45 µm
Sterile Syringe Filter

Millex SLHV033RB [↗](#)

SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

Day 0: Seeding

- 1 18 – 24 hours before transfection, seed 293T cells into a 6 well plate or other format with a density that will make the cells **80 – 95 %** confluent on the day of transfection. Refer to a seeding chart if necessary to seed appropriate density.



Incubate overnight.

Day 1: Transfection

- 3 Change 293T media with fresh DMEM.
- 4 Warm *TransIT*-Lenti Reagent to **Room temperature**.
- 5 Vortex gently before using.
- 6 Gather Opti-Mem, DNA, and packaging mix and refer the table below for the recommended amount of reagents to add based on the format of 293Ts seeded. **Amounts refer to each well of a plate.**



Typically for individual sgRNAs, 2 wells of a 6 well plate per sgRNA will produce enough Lentivirus particles.

For sgRNA libraries (containing up to 50,000 elements), a 15 cm dish can be used. This can be scaled down for smaller libraries.

Culture vessel	48-well plate	24-well plate	12-well plate	6-well plate	10-cm dish	T75 flask	15-cm dish
Surface area	1.0 cm ²	1.9cm ²	3.8 cm ²	9.6 cm ²	59 cm ²	75 cm ²	145cm ²
Complete growth medium	263 µl	0.5 ml	1.0 ml	2.0 ml	10 ml	15 ml	30 ml
Opti-Mem serum-free medium	26 µl	50 µl	100 µl	200 µl	1.0 ml	1.5 ml	3.0 ml
Transfer DNA (1 µg/µl stock)	0.13 µl	0.25 µl	0.5 µl	1.0 µl	5 µl	7.5 µl	15 µl
Packaging DNA Premix (1 µg/µl stock)	0.13 µl	0.25 µl	0.5 µl	1.0 µl	5 µl	7.5 µl	15 µl
<i>TransIT</i> -Lenti Reagent	0.78 µl	1.5 µl	3 µl	6 µl	30 µl	45 µl	90 µl

- 7 Add Opti-MEM into a sterile tube.

8 


In another tube, mix Packaging DNA Premix and DNA.

9 

Add the DNA mix to the Opti-MEM and mix gently.

10 

Add *TransIT*-Lenti Reagent to the mixture and mix gently.

11 

Incubate for  00:10:00 for transfection complexes to form.

12 Add all of the *TransIT*-Lenti:DNA complex mixture to the 293Ts dropwise and gently swirl to mix.

13 

Incubate for 2 days. If a fluorescent marker is included in your DNA, you can check if cells are making virus by checking fluorescence after 24 hours.

Day 3: Harvest

14 With a 12 ml syringe, take up the media/supernatant of the cells.

15 Put a 0.45 μ m filter on the syringe and filter the supernatant into a fresh 15 ml conical tube.



Change the filter if it becomes hard to push. Do not push too hard that bubbles are coming out.


16 Add 1:4 ratio of cold viral precipitation solution (e.g. 0.25 mL viral precipitation solution for 1 mL of viral supernatant).

17 

Mix well by pipetting up and down 10x.

18 

Incubate the viral supernatant at  4 °C for at least  04:00:00 and up to 3 days but no more than 3 days.

- 19 Cool down the centrifuge to 4°C .
- 20  Spin down viral supernatant for 00:30:00 at 1500 x g.
- 21 The pellet will contain the virus.
Resuspend the pellet with 1 ml of your media of choice.
- 22 Virus can be aliquoted and stored at -80°C for long term or 4°C for short term (a few days).



Flash freezing the virus particles in liquid nitrogen may increase the retention of their potency.

Transduction with virus

- 23 Seed iPSC cells so that they will reach 50 % confluency the next day.
- 24 Add virus to cells. The amount to add depends on how concentrated the virus is (adding $\frac{1}{4}$ or $\frac{1}{2}$ of the total produced virus to cells is generally sufficient, see below for typical infection amounts).

step case

sgRNA library (15 cm dish)

A library prepared from a 15 cm dish typically infects 10 million iPSCs in one matrigel-coated T175 flask using 50 % of the produced virus.

- 25 Check next day for fluorescence by microscopy and the next time they are passaged by flow cytometry to check transduction efficiency.
- 26 For sgRNA constructs including puromycin resistance, add 0.8 ug/ml puromycin to select for cells with the sgRNA until they are at least 80 % confluent (typically within 2 passages).

step case

Individual sgRNAs (2 wells of a 6-well plate)

Individual sgRNAs prepared from 2 wells of a 6-well plate typically infect 0.5 million iPSCs in one matrigel-coated well of a 6-well plate using 25 % of the produced virus.

- 25 Check next day for fluorescence by microscopy, and the next time they are passaged by flow cytometry, to check transduction efficiency.
- 26 For sgRNA constructs including puromycin resistance, add 0.8 ug/ml puromycin to select for cells with the sgRNA until they are at least 80 % confluent (typically within 2 passages).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited