



Sep 22, 2019

Transformation using Heat-Shock

[Alba Balletbó](#)¹¹Wageningen University

1

Works for me

[dx.doi.org/10.17504/protocols.io.7j2hkqe](https://doi.org/10.17504/protocols.io.7j2hkqe)

Alba Balletbó



MATERIALS TEXT

Plasmid DNA

SOC Medium

Antibiotic (If applicable)

- 1 Add plasmid DNA (1 µL of miniprep/2 µL of ligation mixture/...µL of cesium-purified DNA).
- 2 Incubate on ice for 30 minutes.
- 3 Heat shock at exactly 42°C for exactly 1:30 minutes, then chill on ice again for 5 minutes.
- 4 Repeat Step 3, once.
- 5 Culture in 1 mL SOC for 1 hour at 37°C and pre-warm the agar plates in the meantime.
- 6 Spin the culture down at 3000rpm for 3 minutes, discard 900 µL and resuspend the pellet in the remaining 100 µL.
- 7 Plate accordingly on agar plates with the appropriate antibiotic, if applicable.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited