

# Q5® Site-Directed Mutagenesis (E0552)

## New England Biolabs

### Abstract

This is the protocol for the Q5® Site-Directed Mutagenesis Kit without competent cells

**Citation:** New England Biolabs Q5® Site-Directed Mutagenesis (E0552). [protocols.io](https://www.protocols.io)

[dx.doi.org/10.17504/protocols.io.chmt45](https://dx.doi.org/10.17504/protocols.io.chmt45)

**Published:** 06 Oct 2014

## Materials

🐛 NEB 5-alpha Competent E.coli (High Efficiency) - 6x0.2 ml [C29871](#) by [New England Biolabs](#)

🐛 Q5 Site-Directed Mutagenesis Kit (Without Competent Cells) - 10 rxns [E0552S](#) by [New England Biolabs](#)

## Protocol

### Exponential Amplification (PCR)

#### Step 1.

Assemble the following reagents in a thin-walled PCR tube.

	25 µl RXN	FINAL CONC.
Q5 Hot Start High-Fidelity 2X Master Mix	12.5 µl	1X
10 µM Forward Primer	1.25 µl	0.5 µM
10 µM Reverse Primer	1.25 µl	0.5 µM
Template DNA (1-25 ng/µl)	1 µl	1-25 ng
Nuclease-free water	9.0 µl	

### ✓ PROTOCOL

#### . [E0552 Q5 PCR Mixture](#)

CONTACT: [New England Biolabs](#)

#### Step 1.1.

Q5 Hot Start High-Fidelity 2X Master Mix **12.5 ul**

#### Step 1.2.

10 µM Forward Primer **1.25ul**

#### Step 1.3.

10 µM Reverse Primer **1.25ul**

#### Step 1.4.

Template DNA (1–25 ng/μl) **1ul**

### Step 1.5.

Nuclease-free water **9ul**

## Exponential Amplification (PCR)

### Step 2.

Mix reagents completely.

## Exponential Amplification (PCR)

### Step 3.

Transfer to a thermocycler and perform the following cycling conditions:

Thermocycling Conditions for a Routine PCR:

STEP	TEMP	TIME
Initial Denaturation	98°C	30 seconds
	98°C	10 seconds
25 Cycles	50–72°C*	10–30 seconds
	72°C	20–30 seconds/kb
Final Extension	72°C	2 minutes
Hold	4–10°C	

## ■ ANNOTATIONS

**New England Biolabs** 09 Oct 2014

\* For a Q5-optimized annealing temperature of mutagenic primers, please use [NEBaseChanger™](#), the online NEB primer design software. For pre-designed, back-to-back primer sets, a  $T_a = T_m + 3$  rule can be applied, but optimization may be necessary.

## Kinase, Ligase & DpnI (KLD) Treatment

### Step 4.

Assemble the following reagents:

## 📄 PROTOCOL

### . [E0552 KLD Mixure](#)

CONTACT: [New England Biolabs](#)

### Step 4.1.

PCR Product **1ul**

### Step 4.2.

2X KLD Reaction Buffer **5ul**

### Step 4.3.

10X KLD Enzyme Mix **1ul**

### Step 4.4.

Nuclease-free Water **3ul**

## Kinase, Ligase & DpnI (KLD) Treatment

### Step 5.

Mix well by pipetting up and down.

## Kinase, Ligase & DpnI (KLD) Treatment

### Step 6.

Incubate at room temperature for 5 minutes.

 **DURATION**

00:05:00

Transformation

### Step 7.

Thaw a **50 µl** aliquot of chemically competent E. coli cells on ice.

■ **ANNOTATIONS**

**New England Biolabs** 26 Jan 2015

NEB 5-alpha Competent E. coli (High Efficiency), [NEB #C2987](#), are recommended

Transformation

### Step 8.

Add **5 µl** of the KLD mix from the "KLD Section" above to the tube of thawed cells.

Transformation

### Step 9.

Carefully flick the tube 4-5 times to mix. **Do not vortex.**

Transformation

### Step 10.

Place the mixture on ice for 30 minutes.

 **DURATION**

00:30:00

Transformation

### Step 11.

Heat shock at 42°C for 30 seconds.

 **DURATION**

00:00:30

Transformation

### Step 12.

Place on ice for 5 minutes.

 **DURATION**

00:05:00

Transformation

### Step 13.

Pipette 950 µl of room temperature SOC into the mixture.

📄 **AMOUNT**

950 µl Additional info:

🧴 **REAGENTS**

 SOC Outgrowth Medium - 100 ml [B9020S](#) by [New England Biolabs](#)

Transformation

### Step 14.

Incubate at 37°C for 60 minutes with shaking (250 rpm).

 **DURATION**

01:00:00

Transformation

### Step 15.

Mix the cells thoroughly by flicking the tube and inverting.

## Transformation

### Step 16.

Spread 50-100 µl onto a selection plate.

## Transformation

### Step 17.

Incubate overnight at 37°C

 **DURATION**

15:00:00

 **NOTES**

**New England Biolabs** 03 Oct 2014

It may be necessary (particularly for simple substitution and deletion experiments) to make a 10- to 40-fold dilution of the transformation mix in SOC prior to plating, to avoid a lawn of colonies