



## GUIDE-seq simplified library preparation protocol (CRISPR/Cas9 off-target cleavage detection)

Nagendra palani<sup>1</sup>

<sup>1</sup>University of Minnesota Genomics Center

[dx.doi.org/10.17504/protocols.io.wikfcccw](https://doi.org/10.17504/protocols.io.wikfcccw)



### ABSTRACT

GUIDE-seq is an experimental method to detect off-target cleavages caused during CRISPR/Cas editing.

See publication <https://doi.org/10.1038/nbt.3117> for details about the method, the protocol used, and the analysis performed.

We identified several areas of improvement in the published protocol for GUIDE-seq NGS library preparation and developed a simplified protocol that is easier to perform & less expensive. The new protocol is also load-and-go compatible with standard sequencing runs on all Illumina instruments (using the kit SBS reagents & default instrument settings) compared to the NBT protocol that requires the libraries to be run on a MiSeq with some tweaks to the software configuration.

### Notable changes:

1. Moved the UMI (Unique Molecular Identifier) from the Y-adapter to be inline with the insert. UMI is now sequenced as part of the sequencing read (R2 on a paired-end run). This eliminates the need to export index reads from the MiSeq.
2. Removed the sample index from the Y-adapter and changed the PCR scheme. Now, only a single adapter is required for ligation to any number of samples, PCR1 is for enrichment of the dsOligo - genome junction, PCR2 is to add sample-specific indexes to the enriched amplicons.

See oligo designs for further information.

3. Replaced sonication based shearing followed by end-repair with a one-pot, one-step enzymatic fragmentation + end-repair kit from New England Biolabs.

### PROTOCOL STATUS

#### Working



We developed this protocol at the University of Minnesota Genomics Center (UMGC) and it is working. UMGc will prepare GUIDE-seq NGS libraries and sequence them as a service.

### GUIDELINES

Time from genomic DNA sample to sequencing-ready library takes 2 days.

Protocol can be performed with 1.5 ml tubes or 96-well plates.

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
TE buffer		
NEBNext Ultra II Q5 Master Mix	E7649	 New England Biolabs
Thermal cycler	View	
NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns	E6177S	 New England Biolabs

### MATERIALS TEXT

Oligonucleotides required:

**TA\_Adaptor\_Top** - /5Phos/CTCACCGCTCTTGTAGSNNNNNNNN CTGTCTCTTATACACATCTCCGAG\*C

**TA\_Adaptor\_Bottom** - CTACAAGAGCGGTGAGT

**dsODN\_Enrich\_Plus** - TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG NNN *GTTTAATTGAGTTGTCATATGTTAATAACG*\*G

**dsODN\_Enrich\_Minus** - TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG NNN *CCGTTATTAACATATGACAACTCAATTAAA*\*C

Replace the italicized sequences with sequences targeting the dsODN that you used in your experiments.

**dsODN\_Enrich\_Adaptor** - GTCTCGTGGGCTCGGAGATGTGTATAAGAGACA\*G

Nextera Indexing Oligo design

**Nextera\_R1** : AATGATACGGCACCACCGAGATCTACAC [i5] TCGTCGGCAGCGT\*C

**Nextera\_R2** : CAAGCAGAAGACGGCATACGAGAT [i7] GTCTCGTGGGCTCG\*G

You will need to replace the [i5] & [i7] with unique index barcodes. You can use the Hamming barcodes from <https://doi.org/10.1371/journal.pone.0036852>. Order the appropriate number of indexing oligos for your samples based on if you want to do combinatorial indexing (cheaper) or unique dual indexing (better with the patterned flowcells).

**Note:** When selecting index barcodes, make sure

- 1) there is no collision with other barcodes that will be run on the same flowcell as your samples
- 2) don't start with or have consecutive G nucleotides if the samples will be run on instruments that use two-color (NextSeq, NovaSeq) or one-color (iSeq) chemistry.

#### SAFETY WARNINGS

NOTE:

The NGS library generated by this protocol is NOT plug-and-play compatible with the guideseq analysis package (<https://github.com/aryeelab/guideseq>).

We have made slight modifications to the guideseq scripts to process these libraries. The output will be the same figure (of off-target sites & their frequencies) as you get with the default package. I will add a comment to this section when the modified package is made available on Github.

#### BEFORE STARTING

Make sure that the genomic DNA is clean (260/230 & 260/280 > 1.8) and 400 ng of DNA is in a volume less than 26 ul. Measure DNA concentration by fluorometry (Picogreen or Qubit).

#### Adaptor formation

1

- Use 0.1x TE to prepare 100 µM stock of all the oligos.
- In a 0.2 ml PCR tube, add the following reagents, vortex well, and spin down. Place on the thermal cycler at **95 °C** for **00:03:00** (heated lid).

**15 µl Oligo TA\_Adaptor\_Top (100 uM)**

**15 µl Oligo TA\_Adaptor\_Bottom (100 uM)**

**70 µl dH2O**

**100 µl Total**

- After incubation, terminate the incubation and let the heat block cool down to room temperature. Remove the tube and place it on the bench for 1 hr.

- The oligos will have annealed and formed the adaptor.

## DNA fragmentation, end-repair, ligation

2

This section of the protocol is from the manual for NEB Ultra II FS kit - Catalog # E7805.

(Refer to <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7805.pdf> : Section for Input > 100 ng for detailed protocol)

- Aliquot **400 ng** of genomic DNA into a well of a 96-well plate and use 1x TE to bring the volume to **26 µl**.
- Add the NEBNext Ultra II FS buffer + enzyme mix to the DNA (follow protocol from the linked manual)
- Set up the following program on the thermal cycler **75 °C heated lid**

**37 °C** for **00:10:00**  
**65 °C** for **00:30:00**  
**4 °C** hold.

### NOTE

This should fragment the library to have a distribution between 300 - 700 bp.

- Add the NEBNext Ultra II ligation mix + enhancer to the well. Add **2.5 µl** of the adaptor to the mixture. Mix well and incubate at **20 °C** for **00:15:00 heated lid off**. Post incubation, add **28.5 µl** of 0.1x TE and bring volume to **100 µl**.
- Perform bead clean-up of the ligation mix using the sample purification beads that came with the kit.  
 Use **20 µl** of beads for 1st bead addition  
**10 µl** of beads for 2nd bead addition (to select for an insert size distribution of 350 - 600 bp).
- Elute in **17 µl** of dH<sub>2</sub>O.

## Enrichment PCR

3

Make 10 uM dilutions of the following primers.

**dsODN\_Enrich\_Plus**  
**dsODN\_Enrich\_Minus**  
**dsODN\_Enrich\_Adaptor**

For each sample set up the following PCRs in a 96-well plate.

*PCR 1.1 - Plus\_Strand\_Enrichment*

**7.5 µl Adaptor ligated DNA**  
**12.5 µl Ultra II Q5 master mix**

2.5 µl dsODN\_Enrich\_Adaptor

2.5 µl dsODN\_Enrich\_Plus

PCR 1.2 - Minus\_Strand\_Enrichment

7.5 µl Adaptor ligated DNA

12.5 µl Ultra II Q5 master mix

2.5 µl dsODN\_Enrich\_Adaptor

2.5 µl dsODN\_Enrich\_Minus

PCR program:

98 °C 00:00:30 Initial denaturation

98 °C 00:00:10 Denaturation

70 °C Gradient -1 °C / cycle 00:00:30 Annealing

72 °C 00:01:00 Extension

7 cycles

98 °C 00:00:10 Denaturation

63 °C 00:00:30 Annealing

72 °C 00:01:00 Extension

13 cycles

4 °C Hold

After end of PCR, use 1 µl of the PCR end-product for Picogreen / Qubit.

Use 2 µl to run the sample on a Tapestation (D5000 high-sensitivity) to check for library size.

Each sample library should be a unimodal distribution between 400 - 1300 bp. All libraries should have an average size distribution within 15% across all samples.

## Indexing PCR

4

Make 10 uM dilutions of the **Nextera\_R1** & **Nextera\_R2** primers.

For each sample, mix

5 ng of PCR 1.1 +

5 ng of PCR 1.2

Set up the following PCR:

10 µl Ultra II Q5 master mix

1 µl Nextera\_R1

1 µl Nextera\_R2

 **10 ng Template**

 **0 µl dH2O** As required

 **20 µl Total**

*PCR Program:*

 **98 °C**  **00:00:30 Initial Denaturation**

 **98 °C**  **00:00:10 Denaturation**

 **65 °C**  **00:00:30 Annealing**

 **72 °C**  **00:01:00 Extension**


10 cycles

 **4 °C Hold**

After end of PCR, use  **1 µl** of the PCR end-product for Picogreen / Qubit.

#### Library pooling, QC, Sequencing

5

- Pool samples by equal mass (assuming the library size distribution for the samples are similar).
- Do a 1.2x SPRI bead clean-up. Elute in  **30 µl 0.1x TE**
- Check concentration of the pool by Picogreen/Qubit and quantify library using Kapa Illumina Library quantification kit. Dilute the pool to appropriate molarity required for the sequencing instrument.
- Load the pool on the sequencer and run with 15% PhiX spike-in for base diversity (works on MiSeq, should work with other Illumina 4-color instruments. Might need more PhiX for NextSeq.)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited