

# Adult Mouse Spleen Dissociation (On ice)

Andrew Potter

## Abstract

Protocol used to dissociate adult (8-10 wk) mouse spleen into single cells. Attained >95% viability, a variety of cell sizes, and ~10 million cells from 12 mg tissue.

**Citation:** Andrew Potter Adult Mouse Spleen Dissociation (On ice). **protocols.io**

[dx.doi.org/10.17504/protocols.io.p2ddqa6](https://doi.org/10.17504/protocols.io.p2ddqa6)

**Published:** 10 May 2018

## Guidelines

### Collagenase Enzyme Mix (two tubes, 1 mL each)

7.5 mg/mL Collagenase A (Sigma, 10103578001)

7.5 mg/mL Collagenase Type 4 (Worthington, CLS-4)

100 µg/mL soybean trypsin inhibitor (Sigma, 10109886001)

125 U DNase (Applchem, A3778)

5 mM CaCl<sub>2</sub>

740 µL DPBS (no Ca, Mg)

---

+12 mg chopped spleen / tube

## Before start

-Set centrifuges to 4° C.

-Make two tubes of 1 mL enzyme mix.

-Make ~25 mL of DPBS/0.04% BSA

## Materials

Red Blood Cell Lysis Buffer Hybri-Max [R7757](#) by [Sigma Aldrich](#)

## Protocol

### Step 1.

Chop tissue coarsely (30 secs) using razor blade on petri dish, on ice.

### Step 2.

Add 12 mg chopped tissue to 1 mL enzyme mix.

### Step 3.

Incubate tube on ice for 10 minutes. Triturate 10X every 2 mins and shake every min.

### Step 4.

After 10 mins of digestion, let tissue chunks settle for 1 min on ice & remove 80% of supernatant with released cells & filter using 70  $\mu$ M filter on 50 mL conical, on ice. Rinse filter with 5 mL ice-cold PBS/0.04% BSA. Leave filter and 50 mL conical on ice, it will be used for the steps as well.

### Step 5.

Add additional 1 mL enzyme mix to tissue chunks.

### Step 6.

Continue to triturate 10x every 2 minutes and shake every minute while incubating on ice, for 10 additional minutes.

### Step 7.

Triturate and add entire volume of cell digestion to 70  $\mu$ M filter on 50 mL conical. Rinse w/5 mL ice-cold PBS/0.04% BSA.

### Step 8.

Transfer flow-through to 15 mL conical. Spin 650 g for five minutes at 4 °C. After spin, remove supernatant (down to 100  $\mu$ L).

### Step 9.

Perform RBC lysis: add 1 mL RBC lysis buffer to cells and triturate 10X. Let sit 3 minutes on ice. Add

10 mL ice-cold PBS/BSA 0.04%

**Step 10.**

Spin 650 g for 5 mins at 4 °C. Remove supernatant and re-suspend in 1 mL PBS/BSA 0.04%.