

# Endothelial paracellular permeability assay

Laura Ruiz Remolina

## Abstract

The human umbilical vein endothelial cell line EA.hy926 (ATCC; Mannasas, VA, USA) was cultured in Dulbecco's Modified Eagle Medium (DMEM) supplemented with 10% fetal calf serum (FCS) in an atmosphere of 5% CO<sub>2</sub>/95% air and 37°C.

EA.hy926 cells (10<sup>5</sup> cells/plate) were seeded in transwell inserts (Corning® Biocoat™ Cell Culture Inserts Collagen, Type I Rat Tail, 24-Well, 3 µm; Corning, NY, USA) and were cultured to confluence at 24 h. The cells were starved overnight and then activated with LPS (1 µg/ml) for 4 h. Upper chamber media, containing LPS (1 µg/ml) and soluble endoglin (500 ng/ml) for their respective treatments, were replaced with FITC-Dextran (40 kDa) at 1 mg/ml in DMEM. The bottom chambers were also replaced with DMEM. After 24 h at 37°C the inserts were removed, and the amount of fluorescence in the bottom chambers was measured using a fluorescence plate reader (Fluoroskan Ascent FL; Thermo Electron Corporation, Waltham, MA, USA).

**Citation:** Laura Ruiz Remolina Endothelial paracellular permeability assay. **protocols.io**

dx.doi.org/10.17504/protocols.io.j8ccrsw

**Published:** 07 Oct 2017

## Guidelines

Avoid repeated freeze-thaw cycles

## Materials

EA.hy926 by [ATCC](#)

Transwell Biocoat Cell culture Insters by [Corning](#)

LPS L3129 by [Sigma-aldrich](#)

✓ Recombinant human endoglin 1097-EN-025 by Contributed by users

DMEM 41966 by [Thermo Fisher Scientific](#)

✓ Fluoroskan Ascent FL by Contributed by users

## Protocol

Culture the cell line EA.hy926 (ATCC; Mannasas, VA, USA) was cultured in DMEM) supplemented with 10% fetal calf serum (FCS) in an atmosphere of 5% CO<sub>2</sub>/95% air and 37°C.

**Step 1.**

Seed EA.hy926 cells (105cells/plate) in transwell inserts (Corning® Biocoat™ Cell Culture Inserts Collagen, Type I Rat Tail, 24-Well, 3 µm; Corning, NY, USA) and cultured it to confluence at 24 h.

**Step 2.**

Starve the cells overnight and activate it with LPS (1 µg/ml) for 4 h.

**Step 3.**

Upper chamber media, containing LPS (1 µg/ml) and soluble endoglin (500 ng/ml) for their respective treatments, were replaced with FITC-Dextran (40 kDa) at 1 mg/ml in DMEM.

**Step 4.**

Replace the bottom chambers of the transwell with DMEM.

**Step 5.**

After 24 h at 37°C remove the inserts and measure the amount of fluorescence in the bottom chambers using a fluorescence plate reader (Fluoroskan Ascent FL; Thermo Electron Corporation, Waltham, MA, USA).

**Step 6.**