

# Algal culture harvest and RNA extraction for RNA-Seq

Alle Lie

## Abstract

**Citation:** Alle Lie Algal culture harvest and RNA extraction for RNA-Seq. **protocols.io**

[dx.doi.org/10.17504/protocols.io.ewjbfcn](https://dx.doi.org/10.17504/protocols.io.ewjbfcn)

**Published:** 26 Apr 2016

## Protocol

### Step 1.

Use sterile techniques to harvest 800 ml of culture, and divide into 4 250 ml centrifuge bottles with 200 ml each.

### Step 2.

Spin the centrifuge bottles at 5000 g for 10 min.

### Step 3.

Carefully decant the supernatant. Use 0.5 ml of RNALater to resuspend the pellet and transfer it to a microcentrifuge tube. Use another 0.25 ml of RNALater to wash the bottle and obtain any remaining cells. Transfer this 0.25 ml into the same microcentrifuge tube.

### Step 4.

Depending on whether you want pseudo-replicates, you can keep the 4 harvested pellet separate or mixed.

Store the harvested cells at or below -20°C until ready for RNA extraction

### Step 5.

Before RNA extraction, the RNALater needs to be removed.

Thaw samples if they are frozen. Centrifuge samples at > 10 k g for 10 min at 4°C. RNALater sometimes affects the buoyancy of the cells, and make them a bit harder to pellet. If you find that the cells are not pelleting, spin harder for a longer time.

Decant supernatant with RNALater

### Step 6.

Extract RNA using a RNeasy Plant Mini Kit (Qiagen) following the manufacturer's instructions

### Step 7.

Remove any DNA from the extracted samples using DNase (Sigma), following manufacturer's instructions

### Step 8.

Clean up the total RNA using a RNA Clean & Concentrator kit (Zymo), following the manufacturer's instructions

Quantify the RNA using any preferred method, e.g. Qubit fluorometer

### Step 9.

Store the RNA at -80°C

It may be ideal to aliquot some RNA out so that you don't need to thaw the entire RNA sample if you need some for any reason