

Nitrite reagent for dissolved cobalt analyses

Randie Bundy

Abstract


This protocol describes the procedures for making a 1.5 mol L⁻¹ nitrite reagent for measuring total dissolved and labile cobalt using cathodic stripping voltammetry.

Citation: Randie Bundy Nitrite reagent for dissolved cobalt analyses. **protocols.io**

dx.doi.org/10.17504/protocols.io.evgbe3w

Published: 17 Apr 2016

Materials

 Sodium nitrite [View](#) by [P212121](#)

✓ MilliQ water by Contributed by users

✓ Chelex-100 resin by Contributed by users

Protocol

Dissolved cobalt analyses

Step 1.

Add 12.9 g of sodium nitrite per 125 mL of Milli-Q, by weighing the sodium nitrite directly into a large trace metal clean bottle.

Dissolved cobalt analyses

Step 2.

Add 125 mL of Milli-Q to the bottle by weight.

Dissolved cobalt analyses

Step 3.

Add clean Chelex-100 resin to the bottle and place it on a shaker overnight to ensure the Chelex mixes throughout the solution.

Dissolved cobalt analyses

Step 4.

The following morning, slowly pour the solution through an HPLC column with a frit in order to filter the chelex and allow the solution to pass into a trace metal clean bottle. Make sure to test the reagent before use to ensure it is clean.