

# E. coli Heat Shock Transformation Version 2

Alex Rajewski

## Abstract

This describes a method to transform a plasmid into homemade DH5 $\alpha$  cells.

**Citation:** Alex Rajewski E. coli Heat Shock Transformation. **protocols.io**

dx.doi.org/10.17504/protocols.io.nwgdgbw

**Published:** 19 Mar 2018

## Materials

- ✓ Dry Block Heater Z3321 by [Promega](#)
- ✓ DH5 $\alpha$  Competent Cells by Contributed by users
- ✓ SOC Media [View](#) by Contributed by users

## Protocol

### Set up

#### Step 1.

Thaw competent cells on ice. This takes about 10 minutes. Also thaw the plasmid to be transformed.

### Set up

#### Step 2.

Warm SOC media to room temperature.

#### AMOUNT

250  $\mu$ l Additional info: SOC media per reaction

#### NOTES

**Alex Rajewski** 19 Mar 2018

LB media can also be used. Anecdotally, we haven't seen a difference between LB and SOC.

### Set up

#### Step 3.

Heat a dry block to 42°C, a shaker-incubator to 37°C, and a cabinet-style incubator to 37°C.

### Set up

#### Step 4.

Label a 1.5mL eppendorf tube for each transformation reaction with the construct name and place it on ice. Label an agar plate with the construct, date, your initials, and plate media type (with antibiotic, if applicable).

## EXPECTED RESULTS

Plate: pYPQ131A 13 March 2018 AR LB+Tet

### Heat Shock

#### Step 5.

Add plasmid to the labeled eppendorf tube on ice.

#### AMOUNT

2 µl Additional info: plasmid

### Heat Shock

#### Step 6.

Add thawed competent cells to the tube with the plasmid and gently swirl with the pipet tip to mix.

#### AMOUNT

45 µl Additional info: competent cells

#### NOTES

**Alex Rajewski** 19 Mar 2018

DO NOT vortex. Competent cells have very specific and fragile cell membranes that can be damaged by vortexing.

### Heat Shock

#### Step 7.

Incubate the mixture on ice for 30 minutes.

#### NOTES

**Alex Rajewski** 19 Mar 2018

This time is fungible. Longer is better, but as short as 5 minutes will also work.

### Heat Shock

#### Step 8.

Place the mixture in the dry block at 42°C for 30 seconds

### Heat Shock

#### Step 9.

Incubate the mixture on ice for 2 minutes.

### Grow Up

## Step 10.

Add SOC to the mixture and shake in an incubator

### AMOUNT

250 µl Additional info: SOC per reaction

Grow Up

## Step 11.

(Optional) Spin the bacterial broth for 30 seconds at 6000 rpm, removed 180µL of supernatant, and gently resuspend the pellet in the remaining media to concentrate.

Plate

## Step 12.

Immerse a plate spreader in 70% ethanol to sterilize and then burn off the excess alcohol.

Plate

## Step 13.

Pipet at least 100µL of the broth onto the labeled agar plate and spread it evenly across the surface of the agar with the sterile spreader. Restilize the spreader.

### NOTES

**Alex Rajewski** 19 Mar 2018

A larger volume can be used, or you can make two plates with different volumes (e.g. 100 and 50µL) to prevent getting colonies that are too close together.

Plate

## Step 14.

Place the prepared plate (lid-down) into a 37°C incubator overnight (16 hours)

### NOTES

**Alex Rajewski** 19 Mar 2018

The plates can also be incubated at room temperature for several days (the weekend).

## Warnings

Competent cells and any materials that touch them should be disposed of in the biohazardous waste container.