



Oct 22, 2019

Total RNA Isolation from Bulk Tissue or Isolated Cells [↗](#)

Brad Godfrey¹¹University of Michigan - Ann Arbor**1** Works for me dx.doi.org/10.17504/protocols.io.8dehs3e**KPMP**Tech. support email: info@kpmp.orgBecky Steck
University of Michigan - Ann Arbor

EXTERNAL LINK

<http://www.kpmp.org>

Procedures

- 1 Transfer biopsy in **1.8 ml** cryotube (polypropylene) in **1 ml** RNAlater
- 2 Use color coding cap inserts to indicate biopsy origin
- 3 Store at **-20 °C** until RNA extraction
- 4 Perform microdissection of glomeruli / tubulointerstitium
- 5 Store dissected tissue portion in less than **50 µl** RNAlater solution and store at **-20 °C** until total RNA extraction
- 6 For total RNA/DNA extraction use "AllPrep DNA/RNA Micro Kit" (# 80284) from Qiagen following respective protocol. Use optional Beta-ME for glomerular compartment.

Note: Elute RNA in **16 µl** H₂O in nucleotide low-binding tubes (Eppendorf #022431021) and store at **-80 °C**.


Note: Keep first flow-through (< 80 bases) for future analysis of small RNAs.



Quality and quantity were measured using Agilent Bioanalyzer 2100

Elute RNA in 16uL

cDNA Generation

- 7 Use NEBNext® Ultra™ II Directional RNA Library Prep Kit for Illumina® (rRNA depletion workflow)1. Use NEBNext® Ultra™ II Directional RNA Library Prep Kit for Illumina® (rRNA depletion workflow)
- 8 1.Store cDNA at  -20 °C

NGS Library Preparation

- 9 Use NovaSeq 6000 to sequence. (The University of Michigan Advanced Genomics Core)1. Use NovaSeq 6000 to sequence. (The University of Michigan Advanced Genomics Core)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited