



May 07,  
2019

Working

## Transformation

Version 4

Rengin Yagmur Akbiyik<sup>1</sup>, Seval Uysalcan<sup>1</sup>

<sup>1</sup>Istanbul Bilgi University

[dx.doi.org/10.17504/protocols.io.2nwgdfc](https://doi.org/10.17504/protocols.io.2nwgdfc)



Rengin Yagmur Akbiyik  
Istanbul Bilgi University



### ABSTRACT

Transformation is a direct modification of the genotype of a cell from a different one by extracellular applications using recombinant DNA techniques. Transformation refers generally to the integration of exogenous DNA into the cell and its integration into the genome. Before starting the transformation process, the gene of interest is ligated with the plasmid vector via the ligase enzyme. This process is called ligation. Generally, nucleic acids or plasmids cannot enter into bacterial cells by themselves. The stimulatory effect is required for that purpose. This means that the cell membranes to be transformed must be pre-arranged. Bacteria that can contain free DNA are called competent bacteria. Some bacteria are highly competitive in normal growth conditions; but some must be treated with chemical or physical methods to gain competitive properties. Competent cells which are cells are generally stimulated with calcium chloride about chemically and the cell membranes are arranged such that the plasmid vector containing the gene of interest can harbor the vector. Some bacterial cells are lasting competent, however; most of them need to be influenced for being competent. There are two ways about this transformation process that;

- Chemical Transformation
- Electrical Transformation (Electroporation)

In general, chemical transformation is used. In chemical way, some special chemicals are used to open the pores which are found in cell membrane. This is because, opening the pores represents the availability of being permeable. Most of the time, divalent ions like Ca<sup>++</sup> are used with assistance of heat-shock.

### STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾	CAS NUMBER ▾	RRID ▾
LB Broth	J106-2KG	Amresco		

Add DNA to competent cells

1

Let stand on ice for 30 min

2

⌚ 00:30:00 On ice

Heat shock

3

⌚ 00:01:00

🔥 42 °C



Grant Bio PHMT Thermoshaker 125-0649



Leave on ice

4 ⌚ 00:01:00

Add LB

5 📄 900 µl for 100 µl cells

📄 800 µl for 200 µl cells



LB Broth

by Amresco

Catalog #: J106-2KG



Shaker -



Put in a shaker

6 🔥 37 °C

⌚ 00:50:00

Centrifuge at 7000 rpm

7 ⌚ 00:01:00



Centrifuge -



Discard LB

8

Resuspend cells and plate them all

9



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited