



Mar 25,  
2019

Working

## Polymerase Chain Reaction - Comparative genomics of *Staphylococcus aureus* associated with subclinical and clinical bovine mastitis (Rocha et al., 2019)

Version 2

Lis Rocha<sup>1</sup>

<sup>1</sup>Universidade Federal de Viçosa

[dx.doi.org/10.17504/protocols.io.zhpf35n](https://doi.org/10.17504/protocols.io.zhpf35n)



Lis Rocha  
Universidade Federal de Viçosa



### ABSTRACT

Comparative genomics of *Staphylococcus aureus* associated with subclinical and clinical bovine mastitis (Rocha et al., 2019)

Many efforts have been made to understand the pathogenesis of bovine mastitis to reduce losses and promote animal welfare. *Staphylococcus aureus* may cause bovine clinical mastitis, but it is mainly associated with subclinical infection, which is usually persistent and can easily reoccur. Here, we conducted a comparative genomic analysis between four strains of *S. aureus* causing subclinical infection (Sau170, 302, 1269, 1364), previously sequenced by our group, and two well-characterized strains causing clinical mastitis (N305 and RF122) to find differences that could be linked to mastitis outcome. A total of 146 virulence-associated genes were compared and no appreciable differences were found between the bacteria. However, several nonsynonymous single nucleotide polymorphisms (SNPs) were identified in genes present in the subclinical strains when compared to RF122, especially in genes encoding host immune evasion and surface proteins. The comparison of orthologous genes using OrthoMCL identified a membrane transporter in the genomes of the bacteria belonging to the subclinical group, but this finding was not confirmed by polymerase chain reaction (PCR) on a collection of field isolates of *S. aureus* associated with clinical or subclinical mastitis. The secreted and surface proteins predicted by different *in silico* tools were compared through multidimensional scaling analysis, revealing a high degree of similarity among the six strains. However, differences were seen in the nucleotide sequences of a gene that codes for a hypothetical protein (cl3309) and a lipoprotein (cl3700). These findings were also analyzed by PCR on DNA extracted from field isolates of *S. aureus*. The lipoprotein, but not the hypothetical protein, was able to separate the clinical isolates from the subclinical ones. These results show that sequence variation among bovine *S. aureus*, and not only the presence/absence of virulence factors, is an important aspect to consider when comparing isolates causing different mastitis outcomes

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS TEXT

50 ng of total DNA, 1U of Taq DNA polymerase Cellco Biotec, 0.2  $\mu$ M of each primer, 0.2 mM deoxynucleotide triphosphate mixture, 1X reaction buffer containing 2.0 mM MgCl<sub>2</sub>, extra 1.0 mM MgCl<sub>2</sub>, and Milli-Q water to increase the reaction volume to a final volume of 25  $\mu$ L.

The extra 1 mM MgCl<sub>2</sub> was excluded from the PCR reaction that contained the primers LipoP-F-CS/LipoP-R-C.

**Table 1 - Primer Sequences for primers used in this Protocol**

cl3309subF	TGTTGTAGGAGGAACAATCC
cl3309subR	TTCTAATGTCAGCAACATGC
cl3309clIF	GCTATTCCTAGATGCACT
cl3309clIR	TTTTAAGTATGACATGAATG
cl3316F	ACGCAAAACCCCTTTACTAGT
cl3316R	GCAACAACCTAGTAGGAGTGA
LipoP-F-CS	GYTTTGCGAAAACGTTAGAYATGTA
LipoP-R-C	TGCCTTCATCATTAAATTGGACCAATC

LipoP-F-CS	GYTTTGCGAAAACGTTAGAYATGTA
LipoP-R-CS	GGTAAAYTCAATGTYCTTATRTCC

#### primers cl3309sub F/R

- 1 Initial denaturation: 95.0 °C for 5 min;
- 2 35 cycles of denaturation at 95.0 °C for 45 s,
- 3 Annealing: 55 °C for 45s
- 4 Extension: 72 °C for 45 s
- 5 final extension at 72.0 °C for 10 min

#### primers cl3316F/R

- 6 initial denaturation: 95.0 °C for 5 min;
- 7 35 cycles of denaturation at 95.0 °C for 45 s,
- 8 Annealing: 55 °C for 45 s
- 9 Extension: 72 °C for 45 s
- 10 final extension at 72.0 °C for 10 min.

#### primers cl3700 - LipoP FCS/RC

- 11 initial denaturation: 95.0 °C for 5 min;
- 12 35 cycles of denaturation at 95.0 °C for 45 s,
- 13 Annealing: 54 °C for 45 s
- 14 Extension: 72 °C for 45 s

15 final extension at 72.0 °C for 10 min.

cl33009cli F/R

16 initial denaturation: 95.0 °C for 5 min;

17 35 cycles of denaturation at 95.0 °C for 45 s,

18 Annealing: 45 °C for 45 s

19 Extension: 72 °C for 30 s

20 final extension at 72.0 °C for 10 min.

primers cl3700 - LipoP FCS/RCS

21 initial denaturation: 95.0 °C for 5 min;

22 35 cycles of denaturation at 95.0 °C for 45 s,

23 Annealing: 50 °C for 45 s

24 Extension: 72 °C for 1min

25 final extension at 72.0 °C for 10 min.

Analyzing the amplified fragments

26 Analyze the amplicons by electrophoresis in 1X Tris-acetate-EDTA on a 1.0% agarose gel and visualize imagen under UV light after staining with 2 mg.ml<sup>-1</sup> ethidium bromide.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited