



Sep 14, 2019

Fluorescence microscope operation procedure

¹¹Northeast Forest University*In Development*

dx.doi.org/10.17504/protocols.io.7cmhiu6

2019 iGEM NEFU_China

Tech. support email: shengyiyenwork@gmail.com

1 Production of temporary fitting

1.1 Add 15ul Bacterial solution to microslide at the clean workbench.

1.2 Add cover slide on it.

2 Use of fluorescence microscopy

2.1 Turn off the light in the room and turn on the microscope mercury lamp.

2.2 Select the appropriate filter according to the luciferin labeled with the sample.

2.3 Place the sample and find the right view. If you want to take a picture, please make sure that the camera is loaded with color film.

2.4 Turn on the selfie device and select the manual file, usually within 0.5 to 10 seconds.

2.5 End of use, turn off all power and record use.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited